

RAB6IP1 siRNA (m): sc-152668

BACKGROUND

RAB6IP1 (rab6-interacting protein 1), also known as DENND5A (DENN domain-containing protein 5A), is a 1,287 amino acid multi-domain protein that belongs to the RAB6IP1 family. Containing a dDENN domain, a DENN domain, a PLAT domain, two RUN domains and a uDENN domain, RAB6IP1 is thought to interact with RAB6A bound to GTP, and may be involved in Rab6-mediated GTPase signaling. The RUN1-PLAT region of Rab6IP1 is required for Rab6-binding. Golgi and endocytic traffic may be linked via Rab6/Rab11 recruitment of Rab6IP1. The RAB6IP1 gene is conserved in chimpanzee, canine, bovine, mouse, rat, zebrafish and *C. elegans*, and maps to human chromosome 11p15.4. With approximately 135 million base pairs and 1,400 genes, chromosome 11 makes up around 4% of human genomic DNA and is considered a gene and disease association dense chromosome. Jervell and Lange-Nielsen syndrome, Jacobsen syndrome, Niemann-Pick disease, hereditary angioedema and Smith-Lemli-Opitz syndrome are all associated with defects in chromosome 11.

REFERENCES

1. Janoueix-Lerosey, I., Jollivet, F., Camonis, J., Marche, P.N. and Goud, B. 1995. Two-hybrid system screen with the small GTP-binding protein Rab6. Identification of a novel mouse GDP dissociation inhibitor isoform and two other potential partners of Rab6. *J. Biol. Chem.* 270: 14801-14808.
2. Grossfeld, P.D., Mattina, T., Lai, Z., Favier, R., Jones, K.L., Cotter, F. and Jones, C. 2004. The 11q terminal deletion disorder: a prospective study of 110 cases. *Am. J. Med. Genet. A* 129A: 51-61.
3. Zehelein, J., Kathoefer, S., Khalil, M., Alter, M., Thomas, D., Brockmeier, K., Ulmer, H.E., Katus, H.A. and Koenen, M. 2006. Skipping of Exon 1 in the KCNQ1 gene causes Jervell and Lange-Nielsen syndrome. *J. Biol. Chem.* 281: 35397-35403.
4. Taylor, T.D., Noguchi, H., Totoki, Y., Toyoda, A., Kuroki, Y., Dewar, K., Lloyd, C., Itoh, T., Takeda, T., Kim, D.W., She, X., Barlow, K.F., Bloom, T., Bruford, E., Chang, J.L., Cuomo, C.A., Eichler, E., Fitzgerald, M.G., Jaffe, D.B., et al. 2006. Human chromosome 11 DNA sequence and analysis including novel gene identification. *Nature* 440: 497-500.
5. Berger, A.C., Salazar, G., Styers, M.L., Newell-Litwa, K.A., Werner, E., Maue, R.A., Corbett, A.H. and Faundez, V. 2007. The subcellular localization of the Niemann-Pick Type C proteins depends on the adaptor complex AP-3. *J. Cell Sci.* 120: 3640-3652.
6. Fernandes, H., Franklin, E., Recacha, R., Houdusse, A., Goud, B. and Khan, A.R. 2009. Structural aspects of Rab6-effector complexes. *Biochem. Soc. Trans.* 37: 1037-1041.
7. Recacha, R., Boulet, A., Jollivet, F., Monier, S., Houdusse, A., Goud, B. and Khan, A.R. 2009. Structural basis for recruitment of Rab6-interacting protein 1 to Golgi via a RUN domain. *Structure* 17: 21-30.

CHROMOSOMAL LOCATION

Genetic locus: Dennd5a (mouse) mapping to 7 F1.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

RAB6IP1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RAB6IP1 shRNA Plasmid (m): sc-152668-SH and RAB6IP1 shRNA (m) Lentiviral Particles: sc-152668-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

RAB6IP1 siRNA (m) is recommended for the inhibition of RAB6IP1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RAB6IP1 gene expression knockdown using RT-PCR Primer: RAB6IP1 (m)-PR: sc-152668-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.