

RABIF siRNA (m): sc-152669

BACKGROUND

The Ras-related superfamily of guanine nucleotide binding proteins includes the R-Ras, Rap, Ral/Rec and Rho/Rab subfamilies, all of which are thought to play an important role in either endocytosis or in biosynthetic protein transport. The process of transporting newly synthesized proteins from the endoplasmic reticulum (ER) to various stacks of the Golgi complex and to secretory vesicles involves the movement of carrier vesicles and requires Rab protein function. Rab proteins are also an integral part of endocytic pathways. RABIF (RAB interacting factor), also known as MSS4, RASGRF3 or RASGRF3, is a 123 amino acid ubiquitously expressed guanine-nucleotide-releasing protein that may participate in vesicular transport. RABIF stimulates GTP-GDP exchange in Sec4 and Rab and binds to a subset of genetically related Rab proteins.

REFERENCES

1. Burton, J.L., Burns, M.E., Gatti, E., Augustine, G.J. and De Camilli, P. 1994. Specific interactions of Mss4 with members of the Rab GTPase subfamily. *EMBO J.* 13: 5547-5558.
2. Miyazaki, A., Sasaki, T., Araki, K., Ueno, N., Imazumi, K., Nagano, F., Takahashi, K. and Takai, Y. 1994. Comparison of kinetic properties between MSS4 and Rab3A GRF GDP/GTP exchange proteins. *FEBS Lett.* 350: 333-336.
3. Yu, H. and Schreiber, S.L. 1995. Cloning, Zn²⁺ binding, and structural characterization of the guanine nucleotide exchange factor human Mss4. *Biochemistry* 34: 9103-9110.
4. Yu, H. and Schreiber, S.L. 1995. Structure of guanine-nucleotide-exchange factor human Mss4 and identification of its Rab-interacting surface. *Nature* 376: 788-791.
5. Müller-Pillasch, F., Zimmerhackl, F., Lacher, U., Schultz, N., Hameister, H., Varga, G., Friess, H., Büchler, M., Adler, G. and Gress, T.M. 1997. Cloning of novel transcripts of the human guanine-nucleotide-exchange factor Mss4: *in situ* chromosomal mapping and expression in pancreatic cancer. *Genomics* 46: 389-396.
6. Thaw, P., Baxter, N.J., Hounslow, A.M., Price, C., Waltho, J.P. and Craven, C.J. 2001. Structure of TCTP reveals unexpected relationship with guanine nucleotide-free chaperones. *Nat. Struct. Biol.* 8: 701-704.
7. Strick, D.J., Francescutti, D.M., Zhao, Y. and Elferink, L.A. 2002. Mammalian suppressor of Sec4 modulates the inhibitory effect of Rab15 during early endocytosis. *J. Biol. Chem.* 277: 32722-32729.
8. Itzen, A., Bleimling, N., Ignatev, A., Pylypenko, O. and Rak, A. 2006. Purification, crystallization and preliminary X-ray crystallographic analysis of mammalian MSS4-Rab8 GTPase protein complex. *Acta Crystallogr. Sect. F Struct. Biol. Cryst. Commun.* 62: 113-116.
9. Zhu, Y.N., Lu, S.M., You, J.F., Zhu, B. and Yu, M.Y. 2009. Novel real-time PCR assay for rapid prenatal diagnosis of Down syndrome: a prospective study of 563 amniocytes. *Clin. Biochem.* 42: 672-675.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Rabif (mouse) mapping to 1 E4.

PRODUCT

RABIF siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RABIF shRNA Plasmid (m): sc-152669-SH and RABIF shRNA (m) Lentiviral Particles: sc-152669-V as alternate gene silencing products.

For independent verification of RABIF (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152669A, sc-152669B and sc-152669C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

RABIF siRNA (m) is recommended for the inhibition of RABIF expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RABIF gene expression knockdown using RT-PCR Primer: RABIF (m)-PR: sc-152669-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.