

RASGEF1C siRNA (m): sc-152709

BACKGROUND

The Ras superfamily of GTPases can be subdivided into the Ras, Rho/Rac, Sar, Rab, ARF and Ran subfamilies and controls multiple aspects of cell function, including cytoskeletal rearrangement, nuclear signaling and cell growth. The Ras superfamily of GTPases function as regulated switches that toggle between a biologically active GTP-bound and an inactive GDP-bound form. This activation is catalyzed by guanine nucleotide exchange factors (GEFs). The RASGEF family of GEFs is highly conserved and the genes encoding the members of this family are present in organisms ranging from nematodes to humans. RASGEF1C (RasGEF domain family member 1C) is a 466 amino acid GEF that contains two RasGEF domains, one at the N-terminus and one at the C-terminus. Due to alternative splicing events, two isoforms exist for RASGEF1C.

REFERENCES

1. Bourne, H.R., Sanders, D.A. and McCormick, F. 1990. The GTPase superfamily: a conserved switch for diverse cell functions. *Nature* 348: 125-132.
2. Boguski, M.S. and McCormick, F. 1993. Proteins regulating Ras and its relatives. *Nature* 366: 643-654.
3. Whitehead, I.P., Khosravi-Far, R., Kirk, H., Trigo-Gonzalez, G., Der, C.J. and Kay, R. 1996. Expression cloning of lsc, a novel oncogene with structural similarities to the Dbl family of guanine nucleotide exchange factors. *J. Biol. Chem.* 271: 18643-18650.
4. Zohn, I.M., Campbell, S.L., Khosravi-Far, R., Rossman, K.L. and Der, C.J. 1998. Rho family proteins and Ras transformation: the RHOad less traveled gets congested. *Oncogene* 17: 1415-1438.
5. Epting, D., Vorwerk, S., Hageman, A. and Meyer, D. 2007. Expression of RASGEF1B in zebrafish. *Gene Expr. Patterns* 7: 389-395.

CHROMOSOMAL LOCATION

Genetic locus: Rasgef1c (mouse) mapping to 11 B1.2.

PRODUCT

RASGEF1C siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RASGEF1C shRNA Plasmid (m): sc-152709-SH and RASGEF1C shRNA (m) Lentiviral Particles: sc-152709-V as alternate gene silencing products.

For independent verification of RASGEF1C (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152709A, sc-152709B and sc-152709C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

RASGEF1C siRNA (m) is recommended for the inhibition of RASGEF1C expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RASGEF1C gene expression knockdown using RT-PCR Primer: RASGEF1C (m)-PR: sc-152709-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.