

▶ RBMS2 siRNA (m): sc-152758

BACKGROUND

RBMS2 is a member of a small family of proteins that bind single stranded DNA or RNA. These proteins are characterized by the presence of two sets of ribonucleoprotein consensus sequence (RNP-CS) that contain conserved motifs, RNP1 and RNP2, and are required for DNA binding. The RBMS proteins have been implicated in such diverse functions as DNA replication, gene transcription, cell cycle progression and apoptosis. RBMS2 (RNA-binding motif, single-stranded-interacting protein 2) is a 407 amino acid protein that contains two RRM (RNA recognition motif) domains and localizes to nucleus. It has been suggested that RBMS2 suppresses Cdc2 kinase and Cdc13 cyclin mutants through the induction of translation of Cdc2. The RBMS2 gene is conserved in chimpanzee, canine, bovine, mouse, rat and zebrafish, and maps to human chromosome 12q13.3.

REFERENCES

1. Skerka, C., Timmann, C., Horstmann, R.D. and Zipfel, P.F. 1992. Two additional human serum proteins structurally related to complement factor H. Evidence for a family of factor H-related genes. *J. Immunol.* 148: 3313-3318.
2. Kanaoka, Y. and Nojima, H. 1994. SCR: novel human suppressors of Cdc2/Cdc13 mutants of *Schizosaccharomyces pombe* harbour motifs for RNA binding proteins. *Nucleic Acids Res.* 22: 2687-2693.
3. Bonaldo, M.F., Lennon, G. and Soares, M.B. 1996. Normalization and subtraction: two approaches to facilitate gene discovery. *Genome Res.* 6: 791-806.
4. Online Mendelian Inheritance in Man, OMIM[™]. 1998. Johns Hopkins University, Baltimore, MD. MIM Number: 602387. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Oppermann, F.S., Gnad, F., Olsen, J.V., Hornberger, R., Greff, Z., Keri, G., Mann, M. and Daub, H. 2009. Large-scale proteomics analysis of the human kinome. *Mol. Cell. Proteomics* 8: 1751-1764.
6. Mayya, V., Lundgren, D.H., Hwang, S.I., Rezaul, K., Wu, L., Eng, J.K., Rodionov, V. and Han, D.K. 2009. Quantitative phosphoproteomic analysis of T cell receptor signaling reveals system-wide modulation of protein-protein interactions. *Sci. Signal.* 2: ra46.

CHROMOSOMAL LOCATION

Genetic locus: Rbms2 (mouse) mapping to 10 D3.

PRODUCT

RBMS2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see RBMS2 shRNA Plasmid (m): sc-152758-SH and RBMS2 shRNA (m) Lentiviral Particles: sc-152758-V as alternate gene silencing products.

For independent verification of RBMS2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-152758A, sc-152758B and sc-152758C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

RBMS2 siRNA (m) is recommended for the inhibition of RBMS2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

RBMS2 (B-4): sc-514918 is recommended as a control antibody for monitoring of RBMS2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor RBMS2 gene expression knockdown using RT-PCR Primer: RBMS2 (m)-PR: sc-152758-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.