

S3-12 siRNA (m): sc-153194

BACKGROUND

S3-12 (adipocyte protein S3-12), also known as PLIN4 (Perilipin-4), is a 1,357 amino acid protein that belongs to the perilipin family. S3-12 contains a repeated 33-amino acid motif also found in adipophilin, and it shares protein sequence identity to both adipophilin and TIP-47 in the COOH terminus, but not to perilipin. S3-12 exhibits sequence similarity to the PAT proteins. Up-regulated during adipocyte differentiation, S3-12 may play a role in triacylglycerol packaging into adipocytes, and functions as a coat protein involved in the biogenesis of lipid storage droplets. S3-12 demonstrates highest expression in skeletal muscle, followed by heart and liver, and lowest expression in adult whole brain. Existing as two alternatively spliced isoforms, the S3-12 gene is conserved in dog, cow and mouse, and maps to human chromosome 19p13.3.

REFERENCES

1. Lu, X., Gruia-Gray, J., Copeland, N.G., Gilbert, D.J., Jenkins, N.A., Londos, C. and Kimmel, A.R. 2001. The murine perilipin gene: the lipid droplet-associated perilipins derive from tissue-specific, mRNA splice variants and define a gene family of ancient origin. *Mamm. Genome* 12: 741-749.
2. Wolins, N.E., Skinner, J.R., Schoenfish, M.J., Tzekov, A., Bensch, K.G. and Bickel, P.E. 2003. Adipocyte protein S3-12 coats nascent lipid droplets. *J. Biol. Chem.* 278: 37713-37721.
3. Dalen, K.T., Schoonjans, K., Ulven, S.M., Weedon-Fekjaer, M.S., Bentzen, T.G., Koutnikova, H., Auwerx, J. and Nebb, H.I. 2004. Adipose tissue expression of the lipid droplet-associating proteins S3-12 and perilipin is controlled by peroxisome proliferator-activated receptor-gamma. *Diabetes* 53: 1243-1252.
4. Grimwood, J., Gordon, L.A., Olsen, A., Terry, A., Schmutz, J., Lamerdin, J., Hellsten, U., Goodstein, D., Couronne, O., Tran-Gyamfi, M., Aerts, A., Altherr, M., Ashworth, L., Bajorek, E., Black, S., Branscomb, E., et al. 2004. The DNA sequence and biology of human chromosome 19. *Nature* 428: 529-535.
5. Kimmel, A.R., Brasaemle, D.L., McAndrews-Hill, M., Sztalryd, C. and Londos, C. 2010. Adoption of PERILIPIN as a unifying nomenclature for the mammalian PAT-family of intracellular lipid storage droplet proteins. *J. Lipid Res.* 51: 468-471.
6. Sone, Y., Yamaguchi, K., Fujiwara, A., Kido, T., Kawahara, K., Ishiwaki, A., Kondo, K., Morita, Y., Tominaga, N. and Otsuka, Y. 2010. Association of lifestyle factors, polymorphisms in adiponectin, perilipin and hormone sensitive lipase, and clinical markers in Japanese males. *J. Nutr. Sci. Vitaminol.* 56: 123-131.
7. Online Mendelian Inheritance in Man, OMIM[™]. 2010. Johns Hopkins University, Baltimore, MD. MIM Number: 613247. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: Plin4 (mouse) mapping to 17 D.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PRODUCT

S3-12 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see S3-12 shRNA Plasmid (m): sc-153194-SH and S3-12 shRNA (m) Lentiviral Particles: sc-153194-V as alternate gene silencing products.

For independent verification of S3-12 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153194A, sc-153194B and sc-153194C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

S3-12 siRNA (m) is recommended for the inhibition of S3-12 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor S3-12 gene expression knockdown using RT-PCR Primer: S3-12 (m)-PR: sc-153194-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.