

SCNM1 siRNA (m): sc-153263

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. SCNM1 (sodium channel modifier 1) is a 230 amino acid protein that contains one matrin-type zinc finger. Localized to the nucleus, SCNM1 is thought to function as an RNA splicing factor that may modify the expression of sodium channel-related proteins. SCNM1 exists as two alternatively spliced isoforms that are encoded by a gene which maps to chromosome 1. Chromosome 1 is the largest human chromosome, spanning about 260 million base pairs and making up 8% of the human genome. Several disorders, including Stickler syndrome, Parkinsons Disease, Gaucher disease, malignant melanoma and Usher syndrome, are caused by defects in genes that localize to chromosome 1.

REFERENCES

1. Sprunger, L.K., Escayg, A., Tallaksen-Greene, S., Albin, R.L. and Meisler, M.H. 1999. Dystonia associated with mutation of the neuronal sodium channel *Scn8a* and identification of the modifier locus SCNM1 on mouse chromosome 3. *Hum. Mol. Genet.* 8: 471-479.
2. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 608095. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. Buchner, D.A., Trudeau, M., George, A.L., Sprunger, L.K. and Meisler, M.H. 2003. High-resolution mapping of the sodium channel modifier SCNM1 on mouse chromosome 3 and identification of a 1.3-kb recombination hot spot. *Genomics* 82: 452-459.
4. Buchner, D.A., Trudeau, M. and Meisler, M.H. 2003. SCNM1, a putative RNA splicing factor that modifies disease severity in mice. *Science* 301: 967-969.
5. Howell, V.M., Jones, J.M., Bergren, S.K., Li, L., Billi, A.C., Avenarius, M.R. and Meisler, M.H. 2007. Evidence for a direct role of the disease modifier SCNM1 in splicing. *Hum. Mol. Genet.* 16: 2506-2516.
6. Howell, V.M., de Haan, G., Bergren, S., Jones, J.M., Culiati, C.T., Michaud, E.J., Frankel, W.N. and Meisler, M.H. 2008. A targeted deleterious allele of the splicing factor SCNM1 in the mouse. *Genetics* 180: 1419-1427.

CHROMOSOMAL LOCATION

Genetic locus: *Scnm1*-ps (mouse) mapping to 3 F2.1.

PRODUCT

SCNM1 siRNA (m) is a pool of 2 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SCNM1 shRNA Plasmid (m): sc-153263-SH and SCNM1 shRNA (m) Lentiviral Particles: sc-153263-V as alternate gene silencing products.

For independent verification of SCNM1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153263A and sc-153263B.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SCNM1 siRNA (m) is recommended for the inhibition of SCNM1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

SCNM1 (A-6): sc-376328 is recommended as a control antibody for monitoring of SCNM1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SCNM1 gene expression knockdown using RT-PCR Primer: SCNM1 (m)-PR: sc-153263-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.