

SerpinA3k siRNA (m): sc-153362

BACKGROUND

The serine proteinase inhibitors (serpins) compose a superfamily of proteins with a diverse set of functions, including the control of blood coagulation, complement activation, programmed cell death and development. Serpins are secreted glycoproteins that contain a stretch of peptide that mimics a true substrate for a corresponding serine protease. SerpinA3k (serine protease inhibitor A3K), also known as contrapsin, SPI-2 or M_{cm}², is a 418 amino acid murine protein that inhibits trypsin-like proteases and belongs to the serpin family. Expressed in liver and secreted in plasma, SerpinA3k is encoded by a gene that maps to murine chromosome 12 E.

REFERENCES

1. Le Cam, A., Pages, G., Auberger, P., Le Cam, G., Leopold, P., Benarous, R. and Glaichenhaus, N. 1987. Study of a growth hormone-regulated protein secreted by rat hepatocytes: cDNA cloning, anti-protease activity and regulation of its synthesis by various hormones. *EMBO J.* 6: 1225-1232.
2. Billingsley, G.D., Walter, M.A., Hammond, G.L. and Cox, D.W. 1993. Physical mapping of four serpin genes: α 1-antitrypsin, α 1-antichymotrypsin, corticosteroid-binding globulin, and protein C inhibitor, within a 280-kb region on chromosome 14q32.1. *Am. J. Hum. Genet.* 52: 343-353.
3. Forsyth, S., Horvath, A. and Coughlin, P. 2003. A review and comparison of the murine α 1-antitrypsin and α 1-antichymotrypsin multigene clusters with the human clade A serpins. *Genomics* 81: 336-345.
4. Horvath, A.J., Forsyth, S.L. and Coughlin, P.B. 2004. Expression patterns of murine antichymotrypsin-like genes reflect evolutionary divergence at the Serpina3 locus. *J. Mol. Evol.* 59: 488-497.
5. Horvath, A.J., Irving, J.A., Rossjohn, J., Law, R.H., Bottomley, S.P., Quinsey, N.S., Pike, R.N., Coughlin, P.B. and Whisstock, J.C. 2005. The murine orthologue of human antichymotrypsin: a structural paradigm for clade A3 serpins. *J. Biol. Chem.* 280: 43168-43178.
6. Winkler, I.G., Hendy, J., Coughlin, P., Horvath, A. and Levesque, J.P. 2005. Serine protease inhibitors serpin1 and serpin3 are down-regulated in bone marrow during hematopoietic progenitor mobilization. *J. Exp. Med.* 201: 1077-1088.

CHROMOSOMAL LOCATION

Genetic locus: Serpina3k (mouse) mapping to 12 E.

PRODUCT

SerpinA3k siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SerpinA3k shRNA Plasmid (m): sc-153362-SH and SerpinA3k shRNA (m) Lentiviral Particles: sc-153362-V as alternate gene silencing products.

For independent verification of SerpinA3k (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153362A, sc-153362B and sc-153362C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SerpinA3k siRNA (m) is recommended for the inhibition of SerpinA3k expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SerpinA3k gene expression knockdown using RT-PCR Primer: SerpinA3k (m)-PR: sc-153362-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.