

SLD5 siRNA (m): sc-153587

BACKGROUND

SLD5 (synthetic lethal with dpb11 mutant-5), also known as GINS4 (GINS complex subunit 4), is a 223 amino acid cytoplasmic and nuclear protein belonging to the GINS4/SLD5 family. Significantly up-regulated in aggressive melanomas, SLD5 is a component of the well-conserved evolutionarily GINS complex, a heterotetramer complex that is also composed of Psf1, Psf2 and Psf3 proteins. The GINS complex plays an important role in the initiation of DNA replication and progression of DNA replication forks. A critical component of the GINS core complex assembly, SLD5 co-localizes with Psf1 and together may cooperate in the proliferation of immature cell populations. Psf1 binds to single-stranded DNA and plays a crucial role in complex function.

REFERENCES

1. Ueno, M., Itoh, M., Kong, L., Sugihara, K., Asano, M. and Takakura, N. 2005. PSF1 is essential for early embryogenesis in mice. *Mol. Cell. Biol.* 25: 10528-10532.
2. Kong, L., Ueno, M., Itoh, M., Yoshioka, K. and Takakura, N. 2006. Identification and characterization of mouse PSF1-binding protein, SLD5. *Biochem. Biophys. Res. Commun.* 339: 1204-1207.
3. Boskovic, J., Coloma, J., Aparicio, T., Zhou, M., Robinson, C.V., Mendez, J. and Montoya, G. 2007. Molecular architecture of the human GINS complex. *EMBO Rep.* 8: 678-684.
4. Kamada, K., Kubota, Y., Arata, T., Shindo, Y. and Hanaoka, F. 2007. Structure of the human GINS complex and its assembly and functional interface in replication initiation. *Nat. Struct. Mol. Biol.* 14: 388-396.
5. Ryu, B., Kim, D.S., Deluca, A.M. and Alani, R.M. 2007. Comprehensive expression profiling of tumor cell lines identifies molecular signatures of melanoma progression. *PLoS ONE* 2: e594.
6. Chang, Y.P., Wang, G., Bermudez, V., Hurwitz, J. and Chen, X.S. 2007. Crystal structure of the GINS complex and functional insights into its role in DNA replication. *Proc. Natl. Acad. Sci. USA* 104: 12685-12690.
7. Nagahama, Y., Ueno, M., Haraguchi, N., Mori, M. and Takakura, N. 2010. PSF3 marks malignant colon cancer and has a role in cancer cell proliferation. *Biochem. Biophys. Res. Commun.* 392: 150-154.
8. MacNeill, S.A. 2010. Structure and function of the GINS complex, a key component of the eukaryotic replisome. *Biochem. J.* 425: 489-500.
9. Pospiech, H., Grosse, F. and Pisani, F.M. 2010. The initiation step of eukaryotic DNA replication. *Subcell. Biochem.* 50: 79-104.

CHROMOSOMAL LOCATION

Genetic locus: Gins4 (mouse) mapping to 8 A2.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

SLD5 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see SLD5 shRNA Plasmid (m): sc-153587-SH and SLD5 shRNA (m) Lentiviral Particles: sc-153587-V as alternate gene silencing products.

For independent verification of SLD5 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as .3 nmol of lyophilized siRNA. These include: sc-153587A, sc-153587B and sc-153587C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

SLD5 siRNA (m) is recommended for the inhibition of SLD5 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor SLD5 gene expression knockdown using RT-PCR Primer: SLD5 (m)-PR: sc-153587-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.