

ST7 siRNA (m): sc-153865

BACKGROUND

ST7 (suppressor of tumorigenicity 7 protein), also known as HELG, RAY1, SEN4, TSG7, ETS7q or FAM4A1, is a 585 amino acid protein localized to the cell membrane. ST7 is ubiquitously expressed, with highest levels found in liver, heart and pancreas. Expressed as seven isoforms produced by alternative splicing events, ST7 may act as a tumor suppressor. The gene that encodes ST7 maps to human chromosome 7, which encodes over 1,000 genes and makes up about 5% of the human genome. Chromosome 7 has been linked to Osteogenesis imperfecta, Pendred syndrome, Lissencephaly, Citrullinemia and Shwachman-Diamond syndrome. The deletion of a portion of the q arm of chromosome 7 is associated with Williams-Beuren syndrome, a condition characterized by mild mental retardation, an unusual comfort and friendliness with strangers and an elfin appearance. Deletions of portions of the q arm of chromosome 7 are also seen in a number of myeloid disorders including cases of acute myelogenous leukemia and myelodysplasia.

REFERENCES

1. Zenklusen, J.C., Conti, C.J. and Green, E.D. 2001. Mutational and functional analyses reveal that ST7 is a highly conserved tumor-suppressor gene on human chromosome 7q31. *Nat. Genet.* 27: 392-398.
2. Brown, V.L., Proby, C.M., Barnes, D.M. and Kelsell, D.P. 2002. Lack of mutations within ST7 gene in tumour-derived cell lines and primary epithelial tumours. *Br. J. Cancer* 87: 208-211.
3. Dong, S.M. and Sidransky, D. 2002. Absence of ST7 gene alterations in human cancer. *Clin. Cancer Res.* 8: 2939-2941.
4. Vincent, J.B., Petek, E., Thevarkunnel, S., Kolozsvari, D., Cheung, J., Patel, M. and Scherer, S.W. 2002. The RAY1/ST7 tumor-suppressor locus on chromosome 7q31 represents a complex multi-transcript system. *Genomics* 80: 283-294.
5. Battle, M.A., Maher, V.M. and McCormick, J.J. 2003. ST7 is a novel low-density lipoprotein receptor-related protein (LRP) with a cytoplasmic tail that interacts with proteins related to signal transduction pathways. *Biochemistry* 42: 7270-7282.
6. Wang, S., Mori, Y., Sato, F., Yin, J., Xu, Y., Zou, T.T., Oлару, A., Kimos, M.C., Perry, K., Selaru, F.M., Deacu, E., Sun, M., Shi, Y.C., Shibata, D., Abraham, J.M., Greenwald, B.D. and Meltzer, S.J. 2003. An LOH and mutational investigation of the ST7 gene locus in human esophageal carcinoma. *Oncogene* 22: 467-470.
7. Sivasundaram, K., Suzuki, H., Seto, M. and Hosokawa, Y. 2003. Mutational analysis of the ST7 gene in human myeloid tumor cell lines. *Oncol. Rep.* 10: 1737-1739.
8. Lu, C., Xu, H.M., Ren, Q., Ao, Y., Wang, Z.N., Ao, X., Jiang, L., Luo, Y. and Zhang, X. 2003. Somatic mutation analysis of p53 and ST7 tumor suppressor genes in gastric carcinoma by DHPLC. *World J. Gastroenterol.* 9: 2662-2665.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: St7 (mouse) mapping to 6 A2.

PRODUCT

ST7 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ST7 shRNA Plasmid (m): sc-153865-SH and ST7 shRNA (m) Lentiviral Particles: sc-153865-V as alternate gene silencing products.

For independent verification of ST7 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153865A, sc-153865B and sc-153865C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ST7 siRNA (m) is recommended for the inhibition of ST7 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ST7 gene expression knockdown using RT-PCR Primer: ST7 (m)-PR: sc-153865-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.