



STAC siRNA (m): sc-153871

BACKGROUND

STAC (SH3 and cysteine rich domain-containing protein), also known as STAC1, is a 402 amino acid protein that contains one SH3 (Src homology 3) domain and one cysteine-rich domain (CRD). Expressed in brain, STAC is a neuron-specific protein that localizes to the cytoplasm and, based on the frequent involvement of SH3 and CRD domains in signal transduction, is believed to play a role in neuron-specific signal transduction. In addition, STAC may be involved in protecting cells from apoptosis. Due to its neuron-specific expression and putative role in signal transduction, STAC may be implicated in a variety of hereditary neurological diseases.

REFERENCES

1. Suzuki, H., Kawai, J., Taga, C., Yaoi, T., Hara, A., Hirose, K., Hayashizaki, Y. and Watanabe, S. 1996. Stac, a novel neuron-specific protein with cysteine-rich and SH3 domains. *Biochem. Biophys. Res. Commun.* 229: 902-909.
2. Kawai, J., Suzuki, H., Hara, A., Hirose, K. and Watanabe, S. 1998. Human and mouse chromosomal mapping of Stac, a neuron-specific protein with an SH3 domain. *Genomics* 47: 140-142.
3. Petek, E., Emberger, W., Kroisel, P.M. and Wagner, K. 1999. Assignment of STAC to human chromosome band 3p22.3 between D3S3718 and D3S1611. *Cytogenet. Cell Genet.* 84: 184-185.
4. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602317. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
5. Hardy, K., Mansfield, L., Mackay, A., Benvenuti, S., Ismail, S., Arora, P., O'Hare, M.J. and Jat, P.S. 2005. Transcriptional networks and cellular senescence in human mammary fibroblasts. *Mol. Biol. Cell* 16: 943-953.
6. Kato, Y., Uzawa, K., Saito, K., Nakashima, D., Kato, M., Nimura, Y., Seki, N. and Tanzawa, H. 2006. Gene expression pattern in oral cancer cervical lymph node metastasis. *Oncol. Rep.* 16: 1009-1014.
7. Satoh, J., Nanri, Y. and Yamamura, T. 2006. Rapid identification of 14-3-3-binding proteins by protein microarray analysis. *J. Neurosci. Methods.* 152: 278-288.

CHROMOSOMAL LOCATION

Genetic locus: Stac (mouse) mapping to 9 F3.

PRODUCT

STAC siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see STAC shRNA Plasmid (m): sc-153871-SH and STAC shRNA (m) Lentiviral Particles: sc-153871-V as alternate gene silencing products.

For independent verification of STAC (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153871A, sc-153871B and sc-153871C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

STAC siRNA (m) is recommended for the inhibition of STAC expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor STAC gene expression knockdown using RT-PCR Primer: STAC (m)-PR: sc-153871-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.