

STK19 siRNA (m): sc-153895

BACKGROUND

The phosphorylation of proteins by protein kinases and protein phosphatases is a key event in most nuclear and cytoplasmic processes. The ability to activate and deactivate proteins via phosphorylation or dephosphorylation is important for cell division, cell differentiation, DNA repair and transcription. STK19 (serine/threonine kinase 19), also known as G11, RP1, D6S60, D6S60E or HLA-RP1, is a 368 amino acid protein that localizes to the nucleus and is a member of the superfamily of serine/threonine protein kinases. Expressed in monocytes, hepatocytes, epithelial cells and T- and B-lymphocytes, STK19 is a protein kinase that can catalytically phosphorylate serine and threonine residues on proteins such as histones and caseins. STK19 functions in an ATP-dependent manner and uses divalent cations, including manganese, as cofactors. Multiple isoforms of STK19 exist due to alternative splicing events.

REFERENCES

1. Shen, L., Wu, L.C., Sanlioglu, S., Chen, R., Mendoza, A.R., Dangel, A.W., Carroll, M.C., Zipf, W.B. and Yu, C.Y. 1994. Structure and genetics of the partially duplicated gene RP located immediately upstream of the complement C4A and the C4B genes in the HLA class III region. Molecular cloning, exon-intron structure, composite retroposon, and breakpoint of gene duplication. *J. Biol. Chem.* 269: 8466-8476.
2. Sargent, C.A., Anderson, M.J., Hsieh, S.L., Kendall, E., Gomez-Escobar, N. and Campbell, R.D. 1994. Characterisation of the novel gene G11 lying adjacent to the complement C4A gene in the human major histocompatibility complex. *Hum. Mol. Genet.* 3: 481-488.
3. Gomez-Escobar, N., Chou, C.F., Lin, W.W., Hsieh, S.L. and Campbell, R.D. 1998. The G11 gene located in the major histocompatibility complex encodes a novel nuclear serine/threonine protein kinase. *J. Biol. Chem.* 273: 30954-30960.
4. Liu, F., Usui, I., Evans, L.G., Austin, D.A., Mellon, P.L., Olefsky, J.M. and Webster, N.J. 2002. Involvement of both G_{q/11} and G_s proteins in gonadotropin-releasing hormone receptor-mediated signaling in L β T2 cells. *J. Biol. Chem.* 277: 32099-32108.
5. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604977. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Wadle, A., Mischo, A., Henrich, P.P., Stenner-Liewen, F., Stenner-Lieven, F., Scherer, C., Imig, J., Petersen, G., Pfreundschuh, M. and Renner, C. 2005. Characterization of Hap/BAG-1 variants as RP1 binding proteins with anti-apoptotic activity. *Int. J. Cancer* 117: 896-904.
7. Roos, C., Dressel, R., Schmidt, B., Günther, E. and Walter, L. 2005. The rat expresses two complement factor C4 proteins, but only one isotype is expressed in the liver. *J. Immunol.* 174: 970-975.

CHROMOSOMAL LOCATION

Genetic locus: Stk19 (mouse) mapping to 17 B1.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

STK19 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see STK19 shRNA Plasmid (m): sc-153895-SH and STK19 shRNA (m) Lentiviral Particles: sc-153895-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

STK19 siRNA (m) is recommended for the inhibition of STK19 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor STK19 gene expression knockdown using RT-PCR Primer: STK19 (m)-PR: sc-153895-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.