

Swiprosin-2 siRNA (m): sc-153965

BACKGROUND

The EF-hand domain is a 12 amino acid loop motif that is commonly found in proteins that participate in calcium-binding events within the cell. EF-hand domains generally exist in a pair that, together, form a stable four-helix bundle that enables the binding of calcium ions. Swiprosin-2, also known as EFHD1 (EF-hand domain-containing protein D1), SWS2, PP3051 or MST133, is a 239 amino acid protein that contains two EF-hand domains and is expressed in a wide variety of tissues, including brain, liver, heart, kidney, testis, ovaries and spleen. Expression of Swiprosin-2, a possible calcium-binding protein, is upregulated during neuronal differentiation, suggesting a role for Swiprosin-2 in brain development and maturation.

REFERENCES

1. Tominaga, M. and Tomooka, Y. 2002. Novel genes cloned from a neuronal cell line newly established from a cerebellum of an adult p53^{-/-} mouse. *Biochem. Biophys. Res. Commun.* 297: 473-479.
2. Lucas, B., Grigo, K., Erdmann, S., Lausen, J., Klein-Hitpass, L. and Ryffel, G.U. 2005. HNF4alpha reduces proliferation of kidney cells and affects genes deregulated in renal cell carcinoma. *Oncogene* 24: 6418-6431.
3. Gifford, J.L., Walsh, M.P. and Vogel, H.J. 2007. Structures and metal-ion-binding properties of the Ca²⁺-binding helix-loop-helix EF-hand motifs. *Biochem. J.* 405: 199-221.
4. Capozzi, F., Luchinat, C., Micheletti, C. and Pontiggia, F. 2007. Essential dynamics of helices provide a functional classification of EF-hand proteins. *J. Proteome Res.* 6: 4245-4255.
5. Online Mendelian Inheritance in Man, OMIM[™]. 2007. Johns Hopkins University, Baltimore, MD. MIM Number: 611617. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Suarez, M.C., Rocha, C.B., Sorenson, M.M., Silva, J.L. and Foguel, D. 2008. Free-energy linkage between folding and calcium binding in EF-hand proteins. *Biophys. J.* 95: 4820-4828.
7. Grigo, K., Wirsing, A., Lucas, B., Klein-Hitpass, L. and Ryffel, G.U. 2008. HNF4 alpha orchestrates a set of 14 genes to down-regulate cell proliferation in kidney cells. *Biol. Chem.* 389: 179-187.

CHROMOSOMAL LOCATION

Genetic locus: Efhd1 (mouse) mapping to 1 D.

PRODUCT

Swiprosin-2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Swiprosin-2 shRNA Plasmid (m): sc-153965-SH and Swiprosin-2 shRNA (m) Lentiviral Particles: sc-153965-V as alternate gene silencing products.

For independent verification of Swiprosin-2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153965A, sc-153965B and sc-153965C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Swiprosin-2 siRNA (m) is recommended for the inhibition of Swiprosin-2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Swiprosin-2 (S-33): sc-100684 is recommended as a control antibody for monitoring of Swiprosin-2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Swiprosin-2 gene expression knockdown using RT-PCR Primer: Swiprosin-2 (m)-PR: sc-153965-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.