



Synaptotagmin X siRNA (m): sc-153976

BACKGROUND

Synaptotagmin X, also known as SytX or synaptotagmin-10 (SYT10), is a 523 amino acid single-pass synaptic vesicle membrane protein that belongs to the synaptotagmin family and contains two C2 domains. Three calcium ions are bound to Synaptotagmin X per subunit using the C2 domains. While it may be involved in calcium-dependent exocytosis of secretory vesicles through calcium and phospholipid binding to the C2 domain, Synaptotagmin X may also serve as calcium sensors in the process of vesicular trafficking and exocytosis. Synaptotagmin X exists as either a homodimer or heterodimer and is only expressed in pancreas, lung and kidney. The gene that encodes Synaptotagmin X contains 64,407 bases and maps to human chromosome 12p11.1. Chromosome 12 is associated with a variety of diseases and afflictions, including hypochondrogenesis, achondrogenesis, Kniest dysplasia, Noonan syndrome and trisomy 12p, which causes facial developmental defects and seizure disorders.

REFERENCES

1. Mikoshiba, K., Fukuda, M., Iwata, K., Kabayama, H. and Mizutani, A. 1999. Role of synaptotagmin, a Ca^{2+} and inositol polyphosphate binding protein, in neurotransmitter release and neurite outgrowth. *Chem. Phys. Lipids* 98: 59-67.
2. Fukuda, M., Kanno, E. and Mikoshiba, K. 1999. Conserved N-terminal cysteine motif is essential for homo- and heterodimer formation of synaptotagmins III, V, VI, and X. *J. Biol. Chem.* 274: 31421-31427.
3. Delgado Carrasco, J., Casanova Morcillo, A., Zabalza Alvillos, M. and Ayala Garces, A. 2001. Achondrogenesis type II-hypochondrogenesis: radiological features. *Case report. An. Esp. Pediatr.* 55: 553-557.
4. Yokoyama, T., Nakatani, S. and Murakami, A. 2003. A case of Kniest dysplasia with retinal detachment and the mutation analysis. *Am. J. Ophthalmol.* 136: 1186-1188.
5. Zhao, E., Li, Y., Fu, X., Zeng, L., Zeng, H., Jin, W., Chen, J., Yin, G., Qian, J., Ying, K., Xie, Y., Zhao, R.C. and Mao, Y. 2003. Cloning and characterization of human synaptotagmin 10 gene. *DNA Seq.* 14: 393-398.
6. Forzano, F., Lituanica, M., Viassolo, A., Superti-Furga, V., Wildhardt, G., Zabel, B. and Faravelli, F. 2007. A familial case of achondrogenesis type II caused by a dominant COL2A1 mutation and "patchy" expression in the mosaic father. *Am. J. Med. Genet. A* 143A: 2815-2820.
7. Wainwright, H. and Beighton, P. 2008. Visceral manifestations of hypochondrogenesis. *Virchows Arch.* 453: 203-207.
8. Lo, F.S., Luo, J.D., Lee, Y.J., Shu, S.G., Kuo, M.T. and Chiou, C.C. 2009. High resolution melting analysis for mutation detection for PTPN11 gene: applications of this method for diagnosis of Noonan syndrome. *Clin. Chim. Acta* 409: 75-77.
9. Benussi, D.G., Costa, P., Zollino, M., Murdolo, M., Petix, V., Carrozzi, M. and Pecile, V. 2009. Trisomy 12p and monosomy 4p: phenotype-genotype correlation. *Genet. Test. Mol. Biomarkers* 13: 199-204.

CHROMOSOMAL LOCATION

Genetic locus: Syt10 (mouse) mapping to 15 E3.

PRODUCT

Synaptotagmin X siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Synaptotagmin X shRNA Plasmid (m): sc-153976-SH and Synaptotagmin X shRNA (m) Lentiviral Particles: sc-153976-V as alternate gene silencing products.

For independent verification of Synaptotagmin X (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-153976A, sc-153976B and sc-153976C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Synaptotagmin X siRNA (m) is recommended for the inhibition of Synaptotagmin X expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Synaptotagmin X gene expression knockdown using RT-PCR Primer: Synaptotagmin X (m)-PR: sc-153976-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.