

TDE1 siRNA (m): sc-154153

BACKGROUND

TDE1, tumor differentially expressed protein 1, is a multi-pass membrane protein also known as serine incorporator 3. Also found on the membranes of the Golgi apparatus within cells, TDE1 is highly expressed in neuronal populations but is also found in thymus, kidney, liver and testis. Expression levels of TDE1 in tumors can be as much as tenfold the amount found in normal tissue of the same type. This increased expression implicates TDE1 as being involved in the cellular transformation from normal to malignant tissue. It is believed TDE1 contributes to oncogenesis by partially protecting cells from serum starvation and etoposide-induced apoptosis. The mechanism through which TDE1 protects cells is poorly understood, but may involve aberrant methylation of TDE1 complexes.

REFERENCES

1. Lebel, M. and Mes-Masson, A.M. 1994. Sequence analysis of a novel cDNA which is overexpressed in testicular tumors from polyomavirus large T-antigen transgenic mice. *DNA Seq.* 5: 31-39.
2. Morrish, D.W., Dakour, J. and Li, H. 1996. Cloning of PL33: a novel probable serpentine membrane receptor associated with human cytotrophoblast and lineage-specific HL-60 cell differentiation. *Mol. Cell. Endocrinol.* 120: 147-151.
3. Bossolasco, M., Lebel, M., Lemieux, N. and Mes-Masson, A.M. 1999. The human TDE gene homologue: localization to 20q13.1-13.3 and variable expression in human tumor cell lines and tissue. *Mol. Carcinog.* 26: 189-200.
4. Nimmrich, I., Erdmann, S., Melchers, U., Finke, U., Hentsch, S., Moyer, M.P., Hoffmann, I. and Müller, O. 2000. Seven genes that are differentially transcribed in colorectal tumor cell lines. *Cancer Lett.* 160: 37-43.
5. Player, A., Gillespie, J., Fujii, T., Fukuoka, J., Dracheva, T., Meerzaman, D., Hong, K.M., Curran, J., Attah, G., Travis, W. and Jen, J. 2003. Identification of TDE2 gene and its expression in non-small cell lung cancer. *Int. J. Cancer* 107: 238-243.
6. Hilton, T.L. and Wang, E.H. 2003. Transcription factor IID recruitment and Sp1 activation. Dual function of TAF1 in cyclin D1 transcription. *J. Biol. Chem.* 278: 12992-13002.
7. Bossolasco, M., Veillette, F., Bertrand, R. and Mes-Masson, A.M. 2006. Human TDE1, a TDE1/TMS family member, inhibits apoptosis *in vitro* and stimulates *in vivo* tumorigenesis. *Oncogene* 25: 4549-4558.
8. Shou, Y., Ma, Z., Lu, T. and Sorrentino, B.P. 2006. Unique risk factors for insertional mutagenesis in a mouse model of XSCID gene therapy. *Proc. Natl. Acad. Sci. USA* 103: 11730-11735.
9. Ariizumi, T., Kawanabe, T., Hatakeyama, K., Sato, S., Kato, T., Tabata, S. and Toriyama, K. 2008. Ultrastructural characterization of exine development of the transient defective exine 1 mutant suggests the existence of a factor involved in constructing reticulate exine architecture from sporopollenin aggregates. *Plant Cell Physiol.* 49: 58-67.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: Serinc3 (mouse) mapping to 2 H3.

PRODUCT

TDE1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TDE1 shRNA Plasmid (m): sc-154153-SH and TDE1 shRNA (m) Lentiviral Particles: sc-154153-V as alternate gene silencing products.

For independent verification of TDE1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154153A, sc-154153B and sc-154153C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TDE1 siRNA (m) is recommended for the inhibition of TDE1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TDE1 gene expression knockdown using RT-PCR Primer: TDE1 (m)-PR: sc-154153-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.