

Tektin 1 siRNA (m): sc-154181

BACKGROUND

The Tektin proteins comprise a family of insoluble, α -helical, filament-forming peptides that interact with Tubulins and, via this interaction, form flagellar and ciliary microtubules. Tektin 1, also known as TEK1, is a 418 amino acid protein that is expressed in testis and contains four evolutionary conserved cysteines and one Tektin-specific peptide sequence. Similar to other members of the Tektin family, Tektin 1 functions as a structural component of microtubules, specifically forming polymers in the walls of microtubules and stabilizing overall microtubule structure. Additionally, Tektin 1 may play a role in spermatogenesis.

REFERENCES

1. Stephens, R.E. and Lemieux, N.A. 1998. Tektins as structural determinants in basal bodies. *Cell Motil. Cytoskeleton* 40: 379-392.
2. Larsson, M., Norrander, J., Gräslund, S., Brundell, E., Linck, R., Stahl, S. and Höög, C. 2000. The spatial and temporal expression of Tekt1, a mouse tektin C homologue, during spermatogenesis suggest that it is involved in the development of the sperm tail basal body and axoneme. *Eur. J. Cell Biol.* 79: 718-725.
3. Xu, M., Zhou, Z., Cheng, C., Zhao, W., Tang, R., Huang, Y., Wang, W., Xu, J., Zeng, L., Xie, Y. and Mao, Y. 2001. Cloning and characterization of a novel human TEK1 gene. *Int. J. Biochem. Cell Biol.* 33: 1172-1182.
4. Iguchi, N., Tanaka, H., Nakamura, Y., Nozaki, M., Fujiwara, T. and Nishimune, Y. 2002. Cloning and characterization of the human Tektin-t gene. *Mol. Hum. Reprod.* 8: 525-530.
5. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 609002. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Clark, A.T., Bodnar, M.S., Fox, M., Rodriguez, R.T., Abeyta, M.J., Firpo, M.T. and Pera, R.A. 2004. Spontaneous differentiation of germ cells from human embryonic stem cells *in vitro*. *Hum. Mol. Genet.* 13: 727-739.
7. Tanaka, H., Iguchi, N., Toyama, Y., Kitamura, K., Takahashi, T., Kaseda, K., Maekawa, M. and Nishimune, Y. 2004. Mice deficient in the axonemal protein Tektin-t exhibit male infertility and immotile-cilium syndrome due to impaired inner arm dynein function. *Mol. Cell. Biol.* 24: 7958-7964.
8. Setter, P.W., Malvey-Dorn, E., Steffen, W., Stephens, R.E. and Linck, R.W. 2006. Tektin interactions and a model for molecular functions. *Exp. Cell Res.* 312: 2880-2896.
9. Amos, L.A. 2008. The tektin family of microtubule-stabilizing proteins. *Genome Biol.* 9: 229.

CHROMOSOMAL LOCATION

Genetic locus: Tekt1 (mouse) mapping to 11 B4.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

Tektin 1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Tektin 1 shRNA Plasmid (m): sc-154181-SH and Tektin 1 shRNA (m) Lentiviral Particles: sc-154181-V as alternate gene silencing products.

For independent verification of Tektin 1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154181A, sc-154181B and sc-154181C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Tektin 1 siRNA (m) is recommended for the inhibition of Tektin 1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Tektin 1 gene expression knockdown using RT-PCR Primer: Tektin 1 (m)-PR: sc-154181-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.