

TIN-Ag siRNA (m): sc-154275

BACKGROUND

TIN-Ag (Tubulointerstitial nephritis antigen), also known as TIN1 or TIN2, is a 476 amino acid secreted protein that contains one SMB (somatomedin-B) domain and localizes to the extracellular matrix. Expressed during development in corneal tissue, kidney cortex and small intestine, TIN-Ag is a basement membrane glycoprotein that interacts with Integrin $\alpha 3$ and Integrin $\alpha 5$ and, via this interaction, mediates the adhesion of proximal tubule epithelial cells. Additionally, TIN-Ag is thought to play a role in the regulation of telomere length, possibly stabilizing the TRF2 (telomeric repeat binding factor 2) complex that is responsible for telomere elongation. Antibodies against TIN-Ag are present in the sera of tubulointerstitial nephritis-affected patients, suggesting a role for TIN-Ag in the pathogenesis of tubulointerstitial nephritis. Two isoforms of TIN-Ag are expressed due to alternative splicing events.

REFERENCES

1. Katz, A., Fish, A.J., Santamaria, P., Nevins, T.E., Kim, Y. and Butkowski, R.J. 1992. Role of antibodies to tubulointerstitial nephritis antigen in human anti-tubular basement membrane nephritis associated with membranous nephropathy. *Am. J. Med.* 93: 691-698.
2. Chen, Y., Krishnamurti, U., Wayner, E.A., Michael, A.F. and Charonis, A.S. 1996. REceptors in proximal tubular epithelial cells for tubulointerstitial nephritis antigen. *Kidney Int.* 49: 153-157.
3. Nelson, T.R., Kim, Y., Michael, A.F., Butkowski, R.J. and Charonis, A.S. 1998. Tubulointerstitial nephritis antigen (TIN-Ag) is expressed in distinct segments of the developing human nephron. *Connect. Tissue Res.* 37: 53-60.
4. Kim, S.H., Kaminker, P. and Campisi, J. 1999. TIN2, a new regulator of telomere length in human cells. *Nat. Genet.* 23: 405-412.
5. Kanwar, Y.S., Kumar, A., Yang, Q., Tian, Y., Wada, J., Kashihara, N. and Wallner, E.I. 1999. Tubulointerstitial nephritis antigen: an extracellular matrix protein that selectively regulates tubulogenesis vs. glomerulogenesis during mammalian renal development. *Proc. Natl. Acad. Sci. USA* 96: 11323-11328.
6. Ikeda, M., Takemura, T., Hino, S. and Yoshioka, K. 2000. Molecular cloning, expression, and chromosomal localization of a human tubulointerstitial nephritis antigen. *Biochem. Biophys. Res. Commun.* 268: 225-230.
7. Zhou, B., Nelson, T.R., Kashtan, C., Gleason, B., Michael, A.F., Vlassi, M. and Charonis, A.S. 2000. Identification of two alternatively spliced forms of human tubulointerstitial nephritis antigen (TIN-Ag). *J. Am. Soc. Nephrol.* 11: 658-668.
8. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 606749. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
9. Ye, J.Z., Donigian, J.R., van Overbeek, M., Loayza, D., Luo, Y., Krutchinsky, A.N., Chait, B.T. and de Lange, T. 2004. TIN2 binds TRF1 and TRF2 simultaneously and stabilizes the TRF2 complex on telomeres. *J. Biol. Chem.* 279: 47264-47271.

CHROMOSOMAL LOCATION

Genetic locus: Tinag (mouse) mapping to 9 D.

PRODUCT

TIN-Ag siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TIN-Ag shRNA Plasmid (m): sc-154275-SH and TIN-Ag shRNA (m) Lentiviral Particles: sc-154275-V as alternate gene silencing products.

For independent verification of TIN-Ag (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154275A, sc-154275B and sc-154275C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TIN-Ag siRNA (m) is recommended for the inhibition of TIN-Ag expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TIN-Ag gene expression knockdown using RT-PCR Primer: TIN-Ag (m)-PR: sc-154275-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.