

TRAM1 siRNA (m): sc-154578

BACKGROUND

TRAM1 (translocation associated membrane protein 1), also known as TRAMP or PNAS-8, is a 374 amino acid protein that localizes to the membrane of the endoplasmic reticulum (ER) and contains one TLC domain. Functioning as a stimulatory protein, TRAM1 influences glycosylation and is required for the translocation of secretory proteins across the ER membrane. The gene encoding TRAM1 maps to human chromosome 8, which consists of nearly 146 million base pairs, houses more than 800 genes and is associated with a variety of diseases and malignancies. Schizophrenia, bipolar disorder, Trisomy 8, Pfeiffer syndrome, congenital hypothyroidism, Waardenburg syndrome and some leukemias and lymphomas are thought to occur as a result of defects in specific genes that maps to chromosome 8.

REFERENCES

1. Görlich, D., Hartmann, E., Prehn, S. and Rapoport, T.A. 1992. A protein of the endoplasmic reticulum involved early in polypeptide translocation. *Nature* 357: 47-52.
2. Walter, P. 1992. Protein translocation. Travelling by TRAM. *Nature* 357: 22-23.
3. Wildenauer, D.B. and Schwab, S.G. 1999. Chromosomes 8 and 10 workshop. *Am. J. Med. Genet.* 88: 239-243.
4. Winter, E. and Ponting, C.P. 2002. TRAM, LAG1 and CLN8: members of a novel family of lipid-sensing domains? *Trends Biochem. Sci.* 27: 381-383.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 605190. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Hegde, R.S. and Kang, S.W. 2008. The concept of translocational regulation. *J. Cell Biol.* 182: 225-232.
7. Kagan, J.C., Su, T., Horng, T., Chow, A., Akira, S. and Medzhitov, R. 2008. TRAM couples endocytosis of Toll-like receptor 4 to the induction of interferon-beta. *Nat. Immunol.* 9: 361-368.

CHROMOSOMAL LOCATION

Genetic locus: Tram1 (mouse) mapping to 1 A3.

PRODUCT

TRAM1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TRAM1 shRNA Plasmid (m): sc-154578-SH and TRAM1 shRNA (m) Lentiviral Particles: sc-154578-V as alternate gene silencing products.

For independent verification of TRAM1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154578A, sc-154578B and sc-154578C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TRAM1 siRNA (m) is recommended for the inhibition of TRAM1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TRAM1 gene expression knockdown using RT-PCR Primer: TRAM1 (m)-PR: sc-154578-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.