

TREML2 siRNA (m): sc-154630

BACKGROUND

TREML2 (triggering receptor expressed on myeloid cells-like 2), also known as C6orf76 or TLT2, is a 321 amino acid single-pass type I membrane protein that contains one Ig-like V-type domain. Expressed in cultured B cells, as well as in T cell leukemia and monocyte leukemia, TREML2 functions as a cell surface receptor that is thought to play a role in the innate and adaptive immune responses. The gene encoding TREML2 maps to a region on human chromosome 6 that is shared by proteins which contain Ig variable (IgV) domains. Chromosome 6 contains 170 million base pairs and comprises nearly 6% of the human genome. Deletion of a portion of the q arm of chromosome 6 is associated with early onset intestinal cancer, suggesting the presence of a cancer susceptibility locus. Additionally, porphyria cutanea tarda, Parkinson's disease, Stickler syndrome and a susceptibility to bipolar disorder are all associated with genes that map to chromosome 6.

REFERENCES

1. Allcock, R.J., Barrow, A.D., Forbes, S., Beck, S. and Trowsdale, J. 2003. The human TREM gene cluster at 6p21.1 encodes both activating and inhibitory single IgV domain receptors and includes NKp44. *Eur. J. Immunol.* 33: 567-577.
2. Online Mendelian Inheritance in Man, OMIM™. 2005. Johns Hopkins University, Baltimore, MD. MIM Number: 609715. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
3. King, R.G., Herrin, B.R. and Justement, L.B. 2006. TREM-like transcript 2 is expressed on cells of the myeloid/granuloid and B lymphoid lineage and is upregulated in response to inflammation. *J. Immunol.* 176: 6012-6021.
4. Wiersinga, W.J., Veer, C.T., Wieland, C.W., Gibot, S., Hooibrink, B., Day, N.P., Peacock, S.J. and van der Poll, T. 2007. Expression profile and function of triggering receptor expressed on myeloid cells-1 during melioidosis. *J. Infect. Dis.* 196: 1707-1716.
5. Hashiguchi, M., Kabori, H., Ritprajak, P., Kamimura, Y., Kozono, H. and Azuma, M. 2008. Triggering receptor expressed on myeloid cell-like transcript 2 (TLT-2) is a counter-receptor for B7-H3 and enhances T cell responses. *Proc. Natl. Acad. Sci. USA* 105: 10495-10500.
6. Molloy, E.J. 2009. Triggering receptor expressed on myeloid cells (TREM) family and the application of its antagonists. *Recent Pat. Antiinfect. Drug Discov.* 4: 51-56.

CHROMOSOMAL LOCATION

Genetic locus: Trem2 (mouse) mapping to 17 C.

PRODUCT

TREML2 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TREML2 shRNA Plasmid (m): sc-154630-SH and TREML2 shRNA (m) Lentiviral Particles: sc-154630-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

TREML2 siRNA (m) is recommended for the inhibition of TREML2 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

TREML2 (G-12): sc-390343 is recommended as a control antibody for monitoring of TREML2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG κ BP-HRP: sc-516102 or m-IgG κ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG κ BP-FITC: sc-516140 or m-IgG κ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TREML2 gene expression knockdown using RT-PCR Primer: TREML2 (m)-PR: sc-154630-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.