

# TRIM67 siRNA (m): sc-154667

## BACKGROUND

TRIM67 (tripartite motif-containing protein 67, TRIM9-like protein) is a 783 amino acid member of the TRIM/RBCC family. Tripartite motif (TRIM) proteins play important roles in a variety of cellular functions including cell proliferation, differentiation, development, oncogenesis, and apoptosis. TRIM gene expression analysis in primary human immune cells seem to suggest the involvement of TRIM proteins in also regulating host antiviral activities. TRIM67 is believed to contain two B box-type zinc fingers, one B30.2/SPRY domain, one COS domain, one fibronectin type-III domain, and one RING-type zinc finger. The gene encoding human TRIM67 is found on chromosome 1q42.2 which is the largest human chromosome spanning about 260 million base pairs and making up 8% of the human genome. There are about 3,000 genes on chromosome 1, and considering the great number of genes there are also a large number of diseases associated with chromosome 1.

## REFERENCES

1. Reymond, A., Meroni, G., Fantozzi, A., Merla, G., Cairo, S., Luzi, L., Riganeli, D., Zanaria, E., Messali, S., Cainarca, S., Guffanti, A., Minucci, S., Pelicci, P.G. and Ballabio, A. 2001. The tripartite motif family identifies cell compartments. *EMBO J.* 20: 2140-2151.
2. Short, K.M. and Cox, T.C. 2006. Subclassification of the RBCC/TRIM superfamily reveals a novel motif necessary for microtubule binding. *J. Biol. Chem.* 281: 8970-8980.
3. Woo, J.S., Suh, H.Y., Park, S.Y. and Oh, B.H. 2006. Structural basis for protein recognition by B30.2/SPRY domains. *Mol. Cell* 24: 967-976.
4. Tao, H., Simmons, B.N., Singireddy, S., Jakkidi, M., Short, K.M., Cox, T.C. and Massiah, M.A. 2008. Structure of the MID1 tandem B-boxes reveals an interaction reminiscent of intermolecular ring heterodimers. *Biochemistry* 47: 2450-2457.
5. Carthagen, L., Bergamaschi, A., Luna, J.M., David, A., Uchil, P.D., Margottin-Goguet, F., Mothes, W., Hazan, U., Transy, C., Pancino, G. and Nisole, S. 2009. Human TRIM gene expression in response to interferons. *PLoS ONE* 4: e4894.
6. Tanji, K., Kamitani, T., Mori, F., Kakita, A., Takahashi, H. and Wakabayashi, K. 2010. TRIM9, a novel brain-specific E3 ubiquitin ligase, is repressed in the brain of Parkinson's disease and dementia with Lewy bodies. *Neurobiol. Dis.* 38: 210-218.
7. Kawai, T. and Akira, S. 2011. Regulation of innate immune signalling pathways by the tripartite motif (TRIM) family proteins. *EMBO Mol. Med.* 3: 513-527.
8. Napolitano, L.M. and Meroni, G. 2012. TRIM family: pleiotropy and diversification through homomultimer and heteromultimer formation. *IUBMB Life* 64: 64-71.

## CHROMOSOMAL LOCATION

Genetic locus: Trim67 (mouse) mapping to 8 E2.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PRODUCT

TRIM67 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see TRIM67 shRNA Plasmid (m): sc-154667-SH and TRIM67 shRNA (m) Lentiviral Particles: sc-154667-V as alternate gene silencing products.

For independent verification of TRIM67 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154667A, sc-154667B and sc-154667C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

TRIM67 siRNA (m) is recommended for the inhibition of TRIM67 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor TRIM67 gene expression knockdown using RT-PCR Primer: TRIM67 (m)-PR: sc-154667-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.