

NET-7 siRNA (m): sc-154725

BACKGROUND

NET-7, also known as TSPAN15 (tetraspanin 15) or TM4SF15 (transmembrane 4 superfamily member 15), is a 294 amino acid multi-pass membrane protein that belongs to the transmembrane 4 superfamily, also known as the tetraspanin family. Members of the tetraspanin family are cell-surface proteins that are characterized by the presence of four hydrophobic domains and mediate signal transduction events that play a role in the regulation of cell development, activation, growth, motility, differentiation, and cancer. Considered molecular facilitators, tetraspanin proteins may regulate vesicle fusion and fission.

REFERENCES

1. Maecker, H.T., Todd, S.C. and Levy, S. 1997. The tetraspanin superfamily: molecular facilitators. *FASEB J.* 11: 428-442.
2. Hemler, M.E. 2003. Tetraspanin proteins mediate cellular penetration, invasion, and fusion events and define a novel type of membrane microdomain. *Annu. Rev. Cell Dev. Biol.* 19: 397-422.
3. Yunta, M. and Lazo, P.A. 2003. Tetraspanin proteins as organisers of membrane microdomains and signalling complexes. *Cell. Signal.* 15: 559-564.
4. Tarrant, J.M., Robb, L., van Spriel, A.B. and Wright, M.D. 2003. Tetraspanins: molecular organisers of the leukocyte surface. *Trends Immunol.* 24: 610-617.
5. Wright, M.D., Moseley, G.W. and van Spriel, A.B. 2004. Tetraspanin microdomains in immune cell signalling and malignant disease. *Tissue Antigens* 64: 533-542.
6. Levy, S. and Shoham, T. 2005. The tetraspanin web modulates immune-signalling complexes. *Nat. Rev. Immunol.* 5: 136-148.
7. Garcia-España, A., Chung, P.J., Sarkar, I.N., Stiner, E., Sun, T.T. and Desalle, R. 2008. Appearance of new tetraspanin genes during vertebrate evolution. *Genomics* 91: 326-334.

CHROMOSOMAL LOCATION

Genetic locus: Tspan15 (mouse) mapping to 10 B4.

PRODUCT

NET-7 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see NET-7 shRNA Plasmid (m): sc-154725-SH and NET-7 shRNA (m) Lentiviral Particles: sc-154725-V as alternate gene silencing products.

For independent verification of NET-7 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154725A, sc-154725B and sc-154725C.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

NET-7 siRNA (m) is recommended for the inhibition of NET-7 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor NET-7 gene expression knockdown using RT-PCR Primer: NET-7 (m)-PR: sc-154725-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.