



# UGP2 siRNA (m): sc-154894

## BACKGROUND

UGP2 (UDP-glucose pyrophosphorylase 2), also known as UPG, UGPP2, UDPGP2 or pHC379, is an evolutionarily conserved protein belonging to the UDPGP type 1 family of proteins. Localizing to the cytoplasm, UGP2 forms homooligomers and is believed to function as a glucosyl donor in cellular metabolic pathways. More specifically, UGP2 catalyzes the transfer of a glucose moiety from glucose-1-phosphate to UTP, producing a diphosphate and UDP-glucose. UDP-glucose is an essential precursor for the synthesis of polysaccharides; in liver and muscle, UDP-glucose is a precursor of glycogen, in liver UDP-glucose is also a precursor of UDP-glucuronate, and in lactating mammary gland UDP-glucose is converted to UDP-galactose and then to lactose.

## REFERENCES

- Shows, T.B., et al. 1978. Assignment of a molecular form of UDP glucose pyrophosphorylase (UGPP2) to chromosome 2 in man. *Cytogenet. Cell Genet.* 22: 215-218.
- Peng, H.L. and Chang, H.Y. 1993. Cloning of a human liver UDP-glucose pyrophosphorylase cDNA by complementation of the bacterial galU mutation. *FEBS Lett.* 329: 153-158.
- Duggleby, R.G., et al. 1996. An improved assay for UDPglucose pyrophosphorylase and other enzymes that have nucleotide products. *Experientia* 52: 568-572.
- Duggleby, R.G., et al. 1996. Sequence differences between human muscle and liver cDNAs for UDPglucose pyrophosphorylase and kinetic properties of the recombinant enzymes expressed in *Escherichia coli*. *Eur. J. Biochem.* 235: 173-179.

## CHROMOSOMAL LOCATION

Genetic locus: Ugp2 (mouse) mapping to 11 A3.1.

## PRODUCT

UGP2 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see UGP2 shRNA Plasmid (m): sc-154894-SH and UGP2 shRNA (m) Lentiviral Particles: sc-154894-V as alternate gene silencing products.

For independent verification of UGP2 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-154894A, sc-154894B and sc-154894C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

UGP2 siRNA (m) is recommended for the inhibition of UGP2 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

UGP2 (B-3): sc-377089 is recommended as a control antibody for monitoring of UGP2 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor UGP2 gene expression knockdown using RT-PCR Primer: UGP2 (m)-PR: sc-154894-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

- Cho, J., et al. 2022. UDP-glucose: a cereblon-dependent glucokinase protein degrader. *Int. J. Mol. Sci.* 23: 9094.
- Cho, J., et al. 2023. Glucokinase variant proteins are resistant to fasting-induced uridine diphosphate glucose-dependent degradation in maturity-onset diabetes of the young type 2 patients. *Int. J. Mol. Sci.* 24: 15842.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.