

VIPL siRNA (m): sc-155109

BACKGROUND

Lectin mannose-binding 1, also designated vesicular integral-membrane protein (VIP36), and lectin mannose-binding 2, also designated ERGIC-53, comprise a family of membrane bound, ubiquitously expressed proteins involved in the selective transport of newly synthesized glycoproteins from the endoplasmic reticulum to the ER-Golgi intermediate compartment. VIPL (VIP36-like protein), also known as LMAN2L (lectin, mannose-binding 2-like), is a 348 amino acid single-pass type I membrane protein that localizes to the endoplasmic reticulum and Golgi apparatus. Containing one L-type lectin-like domain, VIPL is highly expressed in skeletal muscle and kidney, and is found at intermediate levels in heart, liver and placenta, and low levels in brain, thymus, spleen, small intestine and lung. VIPL is suggested to be involved in the regulation of export from the endoplasmic reticulum of a subset of glycoproteins. VIPL may function as a regulator of ERGIC-53. VIPL exists as two alternatively spliced isoforms.

REFERENCES

- Schindler, R., et al. 1993. ERGIC-53, a membrane protein of the ER-Golgi intermediate compartment, carries an ER retention motif. *Eur. J. Cell Biol.* 61: 1-9.
- Kappeler, F., et al. 1994. A dual role for COOH-terminal lysine residues in pre-Golgi retention and endocytosis of ERGIC-53. *J. Biol. Chem.* 269: 6279-6281.
- Hauri, H.P., et al. 2002. Lectins and protein traffic early in the secretory pathway. *Biochem. Soc. Symp.* 69: 73-82.
- Neve, E.P., et al. 2003. VIPL, a VIP36-like membrane protein with a putative function in the export of glycoproteins from the endoplasmic reticulum. *Exp. Cell Res.* 288: 70-83.
- Nufer, O., et al. 2003. Profile-based data base scanning for animal L-type lectins and characterization of VIPL, a novel VIP36-like endoplasmic reticulum protein. *J. Biol. Chem.* 278: 15886-15896.
- Hara-Kuge, S., et al. 2004. The binding of VIP36 and α -amylase in the secretory vesicles via high-mannose type glycans. *Glycobiology* 14: 739-744.
- Kamiya, Y., et al. 2005. Sugar-binding properties of VIP36, an intracellular animal lectin operating as a cargo receptor. *J. Biol. Chem.* 280: 37178-37182.
- Neve, E.P., et al. 2005. Oligomerization and interacellular localization of the glycoprotein receptor ERGIC-53 is independent of disulfide bonds. *J. Mol. Biol.* 354: 556-568.
- Kamiya, Y., et al. 2008. Molecular basis of sugar recognition by the human L-type lectins ERGIC-53, VIPL, and VIP36. *J. Biol. Chem.* 283: 1857-1861.

CHROMOSOMAL LOCATION

Genetic locus: Lman2l (mouse) mapping to 1 B.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

PRODUCT

VIPL siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see VIPL shRNA Plasmid (m): sc-155109-SH and VIPL shRNA (m) Lentiviral Particles: sc-155109-V as alternate gene silencing products.

For independent verification of VIPL (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-155109A, sc-155109B and sc-155109C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

VIPL siRNA (m) is recommended for the inhibition of VIPL expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor VIPL gene expression knockdown using RT-PCR Primer: VIPL (m)-PR: sc-155109-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.