

ZFY1 siRNA (m): sc-155600

BACKGROUND

Zinc-finger proteins contain DNA-binding domains and have a wide variety of functions, most of which encompass some form of transcriptional activation or repression. The majority of zinc-finger proteins contain a Krüppel-type DNA binding domain and a KRAB domain, which is thought to interact with KAP1, thereby recruiting histone modifying proteins. ZFY1 (zinc finger protein 1, Y linked), also known as ZFY1 is a 782 amino acid nuclear protein belonging to the Krüppel C₂H₂-type zinc-finger protein family and the ZFX/ZFY subfamily. Containing 13 C₂H₂-type zinc fingers, ZFY1 is expressed in the genital ridge and adult testis and may be a probable transcription activator. The gene encoding ZFY1 maps to mouse chromosome Y A1.

REFERENCES

1. Ashworth, A., Swift, S. and Affara, N. 1989. Sequence of cDNA for murine Zfy-1, a candidate for Tdy. *Nucleic Acids Res.* 17: 2864.
2. Hayashida, H., Kuma, K. and Miyata, T. 1992. Interchromosomal gene conversion as a possible mechanism for explaining divergence patterns of ZFY-related genes. *J. Mol. Evol.* 35: 181-183.
3. Pamilo, P. and Bianchi, N.O. 1993. Evolution of the Zfx and Zfy genes: rates and interdependence between the genes. *Mol. Biol. Evol.* 10: 271-281.
4. Zambrowicz, B.P., Zimmermann, J.W., Harendza, C.J., Simpson, E.M., Page, D.C., Brinster, R.L. and Palmiter, R.D. 1994. Expression of a mouse Zfy-1/lacZ transgene in the somatic cells of the embryonic gonad and germ cells of the adult testis. *Development* 120: 1549-1559.
5. Zambrowicz, B.P., Findley, S.D., Simpson, E.M., Page, D.C. and Palmiter, R.D. 1994. Characterization of the murine Zfy1 and Zfy2 promoters. *Genomics* 24: 406-408.
6. Mahaffey, C.L., Bayleran, J.K., Yeh, G.Y., Lee, T.C., Page, D.C. and Simpson, E.M. 1997. Intron/exon structure confirms that mouse Zfy1 and Zfy2 are members of the ZFY gene family. *Genomics* 41: 123-127.
7. Nagamine, C.M. and Carlisle, C. 1997. Zfy2/1 fusion gene fails to replicate Zfy1 expression pattern in fetal gonads. *Genomics* 43: 397-398.
8. Johnston, C.M., Shimeld, S.M. and Sharpe, P.T. 1998. Molecular evolution of the ZFY and ZNF6 gene families. *Mol. Biol. Evol.* 15: 129-137.
9. Hansen, M.A., Nielsen, J.E., Tanaka, M., Almstrup, K., Skakkebaek, N.E. and Leffers, H. 2006. Identification and expression profiling of 10 novel spermatid expressed CYPT genes. *Mol. Reprod. Dev.* 73: 568-579.

CHROMOSOMAL LOCATION

Genetic locus: Zfy1 (mouse) mapping to Y A1.

PRODUCT

ZFY1 siRNA (m) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ZFY1 shRNA Plasmid (m): sc-155600-SH and ZFY1 shRNA (m) Lentiviral Particles: sc-155600-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ZFY1 siRNA (m) is recommended for the inhibition of ZFY1 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ZFY1 gene expression knockdown using RT-PCR Primer: ZFY1 (m)-PR: sc-155600-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.