SANTA CRUZ BIOTECHNOLOGY, INC.

CA VIII (E-4): sc-166626



BACKGROUND

Carbonic anhydrases (CAs) are members of a large family of zinc metalloenzymes responsible for catalyzing the reversible hydration of carbon dioxide. CAs show extensive diversity in their distribution and subcellular localization. They are involved in a variety of biological processes, including calcification, bone resorption, respiration, acid-base balance and the formation of aqueous humor, saliva, gastric juice and cerebrospinal fluid. CA VIII, also referred to as carbonic anhydrase-related protein VIII (CA-RP VIII), is a member of the carbonic anhydrase family that lacks the Zn-binding motif essential for carbonic anhydrase activity. For this reason, CA VIII does not exhibit catalytic activity but instead may be important in synaptic vesicle formation and transport. In addition, CA VIII may be involved in the invasiveness of nonsmall cell lung carcinomas and may also play a role in the growth of colon cancer cells.

REFERENCES

- 1. Bataller, L., et al. 2004. Carbonic anhydrase-related protein VIII: autoantigen in paraneoplastic cerebellar degeneration. Ann. Neurol. 56: 575-579.
- Morimoto, K., et al. 2005. Overexpression of carbonic anhydrase-related protein XI promotes proliferation and invasion of gastrointestinal stromal tumors. Virchows Arch. 447: 66-73.
- Jiao, Y., et al. 2005. Carbonic anhydrase-related protein VIII deficiency is associated with a distinctive lifelong gait disorder in waddles mice. Genetics 171: 1239-1246.

CHROMOSOMAL LOCATION

Genetic locus: CA8 (human) mapping to 8q12.1; Car8 (mouse) mapping to 4 A1.

SOURCE

CA VIII (E-4) is a mouse monoclonal antibody raised against amino acids 1-100 mapping at the N-terminus of CA VIII of human origin.

PRODUCT

Each vial contains 200 μg IgG_{2a} kappa light chain in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

CA VIII (E-4) is available conjugated to agarose (sc-166626 AC), 500 μ g/ 0.25 ml agarose in 1 ml, for IP; to HRP (sc-166626 HRP), 200 μ g/ml, for WB, IHC(P) and ELISA; to either phycoerythrin (sc-166626 PE), fluorescein (sc-166626 AF546; Alexa Fluor[®] 488 (sc-166626 AF488), Alexa Fluor[®] 546 (sc-166626 AF546), Alexa Fluor[®] 594 (sc-166626 AF594) or Alexa Fluor[®] 647 (sc-166626 AF647), 200 μ g/ml, for WB (RGB), IF, IHC(P) and FCM; and to either Alexa Fluor[®] 680 (sc-166626 AF680) or Alexa Fluor[®] 790 (sc-166626 AF790), 200 μ g/ml, for Near-Infrared (NIR) WB, IF and FCM.

Alexa Fluor® is a trademark of Molecular Probes, Inc., Oregon, USA

STORAGE

Store at 4° C, **D0 NOT FREEZE**. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

APPLICATIONS

CA VIII (E-4) is recommended for detection of CA VIII of mouse, rat and human origin by Western Blotting (starting dilution 1:100, dilution range 1:100-1:1000), immunoprecipitation [1-2 μ g per 100-500 μ g of total protein (1 ml of cell lysate)], immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500), immunohistochemistry (including paraffin-embedded sections) (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

Suitable for use as control antibody for CA VIII siRNA (h): sc-62038, CA VIII siRNA (m): sc-62039, CA VIII shRNA Plasmid (h): sc-62038-SH, CA VIII shRNA Plasmid (m): sc-62039-SH, CA VIII shRNA (h) Lentiviral Particles: sc-62038-V and CA VIII shRNA (m) Lentiviral Particles: sc-62039-V.

Molecular Weight of CA VIII: 33 kDa.

Positive Controls: C32 whole cell lysate: sc-2205, SK-MEL-28 cell lysate: sc-2236 or mouse brain extract: sc-2253.

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunoprecipitation: use Protein A/G PLUS-Agarose: sc-2003 (0.5 ml agarose/2.0 ml). 3) Immunofluorescence: use m-IgGκ BP-FITC: sc-516140 or m-IgGκ BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz[®] Mounting Medium: sc-24941 or UltraCruz[®] Hard-set Mounting Medium: sc-359850. 4) Immunohistochemistry: use m-IgGκ BP-HRP: sc-516102 with DAB, 50X: sc-24982 and Immunohistomount: sc-45086, or Organo/Limonene Mount: sc-45087.

DATA





CA VIII (E-4): sc-166626. Western blot analysis of CA VIII expression in C32 (${\rm A})$ and SK-MEL-28 (${\rm B})$ whole cell lysates.

CA VIII (E-4): sc-166626. Immunoperoxidase staining of formalin fixed, paraffin-embedded human (A) and mouse (B) cerebellum tissue showing cytoplasmic staining of Purkinje cells and neuropil staining in molecular layer.

SELECT PRODUCT CITATIONS

1. White, J.J., et al. 2016. An optimized surgical approach for obtaining stable extracellular single-unit recordings from the cerebellum of head-fixed behaving mice. J. Neurosci. Methods 262: 21-31.

RESEARCH USE

For research use only, not for use in diagnostic procedures.