

G_{αz} (h4): 293T Lysate: sc-170076

BACKGROUND

Heterotrimeric G proteins function to relay information from cell surface receptors to intracellular effectors. Each of a very broad range of receptors specifically detects an extracellular stimulus (a photon, pheromone, odorant, hormone or neurotransmitter) while the effectors (i.e., adenylyl cyclase), which act to generate one or more intracellular messengers, are less numerous. In mammals, G protein α , β and γ polypeptides are encoded by at least 16, 4 and 7 genes, respectively. Most interest in G proteins has been focused on their α subunits, since these proteins bind and hydrolyze GTP and most obviously regulate the activity of the best studied effectors. Four distinct classes of G $_{\alpha}$ subunits have been identified; these include G $_s$, G $_i$, G $_q$ and G $_{\alpha 12/13}$. The G $_i$ class comprises all the known α subunits that are susceptible to pertussis toxin modifications, including G $_{\alpha i-1}$, G $_{\alpha i-2}$, G $_{\alpha i-3}$, G $_{\alpha o}$, G $_{\alpha 11}$, G $_{\alpha 12}$, G $_{\alpha z}$ and G $_{\alpha gust}$. Of these, the three G $_{\alpha i}$ subtypes function to open atrial potassium channels.

REFERENCES

1. Jones, D.T., Masters, S.B., Bourne, H.R. and Reed, R.R. 1990. Biochemical characterization of three stimulatory GTP-binding proteins. The large and small forms of G $_s$ and the olfactory-specific G-protein, G $_{olf}$. J. Biol. Chem. 265: 2671-2676.
2. Simon, M.I., Strathmann, M.P. and Gautam, N. 1991. Diversity of G proteins in signal transduction. Science 252: 802-808.
3. Cali, J.J., Balcueva, E.A., Rybalkin, I. and Robishaw, J.D. 1992. Selective tissue distribution of G protein γ subunits, including a new form of the γ subunits identified by cDNA cloning. J. Biol. Chem. 267: 24023-24027.
4. McLaughlin, S.K., McKinnon, P.J. and Margolskee, R.F. 1992. Gustducin is a taste cell-specific G protein closely related to the transducins. Nature 357: 563-569.
5. von Weizsäcker, E., Strathman, M.P. and Simon, M.I. 1992. Diversity among the β subunits of heterotrimeric GTP-binding proteins: characterization of a novel β -subunit cDNA. Biochem. Biophys. Res. Commun. 183: 350-356.
6. Conklin, B.R. and Bourne, H.R. 1993. Structural elements of G $_{\alpha}$ subunits that interact with G $_{\beta\gamma}$, receptors, and effectors. Cell 73: 631-641.

CHROMOSOMAL LOCATION

Genetic locus: GNAZ (human) mapping to 22q11.22.

PRODUCT

G $_{\alpha z}$ (h): 293T Lysate represents a lysate of human G $_{\alpha z}$ transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

G $_{\alpha z}$ (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive G $_{\alpha z}$ antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.