

CNPase (h2): 293 Lysate: sc-170132

BACKGROUND

2',3'-cyclic nucleotide 3'-phosphodiesterase (CNPase) is a membrane-bound enzyme that can link Tubulin to membranes and may regulate cytoplasmic microtubule distribution. CNPase acts as a microtubule-associated protein by promoting microtubule assembly; this activity resides in the C-terminus of the enzyme. CNPase is firmly associated with Tubulin from brain tissue and thyroid cells and can be found at high concentrations in central nervous system Myelin and in the outer segments of photoreceptors in the retina. Acute lead intoxication leads to disturbances in CNPase activity and the morphology of Myelin.

REFERENCES

1. Sprinkle, T.J., et al. 1987. Monoclonal antibody production to human and bovine 2',3'-cyclic nucleotide 3'-phosphodiesterase (CNPase): high-specificity recognition in whole brain acetone powders and conservation of sequence between CNP1 and CNP2. *Brain Res.* 426: 349-357.
2. Vogel, U., et al. 1988. Molecular structure, localization, and possible functions of the Myelin-associated enzyme 2',3'-cyclic nucleotide 3'-phosphodiesterase. *J. Neurochem.* 50: 1667-1677.
3. Dabrowska-Bouta, B., et al. 2000. Acute lead intoxication *in vivo* affects Myelin membrane morphology and CNPase activity. *Exp. Toxicol. Pathol.* 52: 257-263.
4. Bifulco, M., et al. 2002. 2',3'-cyclic nucleotide 3'-phosphodiesterase: a membrane-bound, microtubule-associated protein and membrane anchor for Tubulin. *Proc. Natl. Acad. Sci. USA* 99: 1807-1812.
5. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 123830. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: CNP (human) mapping to 17q21.2.

PRODUCT

CNPase (h2): 293 Lysate represents a lysate of human CNPase transfected 293 cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

APPLICATIONS

CNPase (h2): 293 Lysate is suitable as a Western Blotting positive control for human reactive CNPase antibodies. Recommended use: 10-20 µl per lane.


Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

CNPase (B-9): sc-166254 is recommended as a positive control antibody for Western Blot analysis of enhanced human CNPase expression in CNPase transfected 293 cells (starting dilution 1:100, dilution range 1:100-1:1,000).


RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
 1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



CNPase (B-9): sc-166254. Western blot analysis of CNPase expression in non-transfected: sc-110760 (A) and human CNPase transfected: sc-170132 (B) 293 whole cell lysates.



CNPase (G-6): sc-166063. Western blot analysis of CNPase expression in non-transfected: sc-110760 (A) and human CNPase transfected: sc-170132 (B) 293 whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.