

β Enolase (h2): 293T Lysate: sc-170445

BACKGROUND

Enolases have been characterized as highly conserved cytoplasmic glycolytic enzymes that may be involved in differentiation. Three isoenzymes have been identified: α Enolase, β Enolase and γ Enolase. α Enolase expression has been detected in most tissues, whereas β Enolase is expressed predominantly in muscle tissue and γ Enolase is detected only in nervous tissue. These isoforms exist as both homodimers and heterodimers, and they play a role in converting phosphoglyceric acid to phosphoenolpyruvic acid in the glycolytic pathway. β Enolase, also known as Enolase 3 or MSE (muscle-specific Enolase), localizes to the cytoplasm and is expressed as a homodimer or a heterodimer with α Enolase in adult skeletal muscle. Mutations in the gene encoding β Enolase may result in glycogenesis type XIII (muscle-specific β Enolase deficiency), a disorder characterized by fatigability, muscle weakness and exercise-induced myalgia (or muscle pain).

REFERENCES

1. Chen, S.H. and Giblett, E.R. 1976. Enolase: human tissue distribution and evidence for three different loci. *Ann. Hum. Genet.* 39: 277-280.
2. Whitehead, M.C., et al. 1982. Synapse formation is related to the onset of neuron-specific enolase immunoreactivity in the avian auditory and vestibular systems. *Dev. Neurosci.* 5: 298-307.
3. Verma, M., et al. 1994. DNA sequences encoding Enolase are remarkably conserved from yeast to mammals. *Life Sci.* 55: 893-899.
4. Keller, A., et al. 1994. Coexpression of α and γ Enolase genes in neurons of adult rat brain. *J. Neurosci. Res.* 38: 493-504.
5. Zhang, E., et al. 1997. Mechanism of Enolase: the crystal structure of asymmetric dimer Enolase-2-phospho-D glycerate/Enolase-phosphoenolpyruvate at 2.0 Å resolution. *Biochemistry* 36: 12526-12534.
6. Deloulme, J.C., et al. 1997. A comparative study of the distribution of α and γ Enolase subunits in cultured rat neural cells and fibroblasts. *Int. J. Dev. Neurosci.* 15: 183-194.
7. Sensenbrenner, M., et al. 1997. Expression of two neuronal markers, growth-associated protein 43 and neuron-specific Enolase, in rat glial cells. *J. Mol. Med.* 75: 653-663.
8. Comi, G.P., et al. 2001. β-Enolase deficiency, a new metabolic myopathy of distal glycolysis. *Ann. Neurol.* 50: 202-207.
9. Online Mendelian Inheritance in Man, OMIM[™]. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 131370. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>

CHROMOSOMAL LOCATION

Genetic locus: ENO3 (human) mapping to 17p13.2.

PRODUCT

β Enolase (h2): 293T Lysate represents a lysate of human β Enolase transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

β Enolase (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive β Enolase antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

β Enolase (XX-10): sc-100811 is recommended as a positive control antibody for Western Blot analysis of enhanced human β Enolase expression in β Enolase transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker[™] Molecular Weight Standards: sc-2035, UltraCruz[®] Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



β Enolase (XX-10): sc-100811. Western blot analysis of β Enolase expression in non-transfected: sc-117752 (A) and human β Enolase transfected: sc-170445 (B) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.