

CRTAP (h2): 293T Lysate: sc-170519

BACKGROUND

CRTAP (cartilage associated protein), also known as CASP or LEPREL3 (leprecan-like 3), is a secreted protein localizing to the extracellular space that plays a role in collagen post-translational modifications, extracellular fibril assembly and intracellular trafficking. CRTAP is widely expressed with predominant expression in articular chondrocytes. It contains a signal peptide and a tetra-ricopeptide-like helical domain and is essential for normal bone formation. In the endoplasmic reticulum (ER), CRTAP forms a complex with Gros1 and CyPB (cyclophilin B) and is required for the efficient 3-hydroxylation of target prolyl residues in collagen type I molecules, the major structural proteins of skin and bone. Mutations in the gene encoding CRTAP can lead to autosomal recessive osteogenesis imperfecta (OI) type 7 and type 2B. OI, also known as brittle bone disease, is characterized by bone fragility and susceptibility to fractures. OI type 7 is a mild form of this disorder, while OI type 2B is a neonatal lethal condition.

REFERENCES

1. Castagnola, P., Gennari, M., Morello, R., Tonachini, L., Marin, O., Gaggero, A. and Cancedda, R. 1997. Cartilage associated protein (CASP) is a novel developmentally regulated chick embryo protein. *J. Cell Sci.* 110: 1351-1359.
2. Morello, R., Tonachini, L., Monticone, M., Viggiano, L., Rocchi, M., Cancedda, R. and Castagnola, P. 1999. cDNA cloning, characterization and chromosome mapping of CRTAP encoding the mouse cartilage associated protein. *Matrix Biol.* 18: 319-324.
3. Tonachini, L., Morello, R., Monticone, M., Skaug, J., Scherer, S.W., Cancedda, R. and Castagnola, P. 1999. cDNA cloning, characterization and chromosome mapping of the gene encoding human cartilage associated protein (CRTAP). *Cytogenet. Cell Genet.* 87: 191-194.
4. Barnes, A.M., Chang, W., Morello, R., Cabral, W.A., Weis, M., Eyre, D.R., Leikin, S., Makareeva, E., Kuznetsova, N., Uveges, T.E., Ashok, A., Flor, A.W., Mulvihill, J.J., Wilson, P.L., Sundaram, U.T., Lee, B. and Marini, J.C. 2006. Deficiency of cartilage-associated protein in recessive lethal osteogenesis imperfecta. *N. Engl. J. Med.* 355: 2757-2764.
5. Morello, R., Bertin, T.K., Chen, Y., Hicks, J., Tonachini, L., Monticone, M., Castagnola, P., Rauch, F., Glorieux, F.H., Vranka, J., Bächinger, H.P., Pace, J.M., Schwarze, U., Byers, P.H., Weis, M., Fernandes, R.J., Eyre, D.R., Yao, Z., Boyce, B.F. and Lee, B. 2006. CRTAP is required for prolyl 3-hydroxylation and mutations cause recessive osteogenesis imperfecta. *Cell* 127: 291-304.
6. Martin, E. and Shapiro, J.R. 2007. Osteogenesis imperfecta: epidemiology and pathophysiology. *Curr. Osteoporos. Rep.* 5: 91-97.
7. Kwan, T., Benovoy, D., Dias, C., Gurd, S., Serre, D., Zuzan, H., Clark, T.A., Schweitzer, A., Staples, M.K., Wang, H., Blume, J.E., Hudson, T.J., Sladek, R. and Majewski, J. 2007. Heritability of alternative splicing in the human genome. *Genome Res.* 17: 1210-1218.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

CHROMOSOMAL LOCATION

Genetic locus: CRTAP (human) mapping to 3p22.3.

PRODUCT

CRTAP (h2): 293T Lysate represents a lysate of human CRTAP transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

CRTAP (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive CRTAP antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

CRTAP (E-1): sc-393136 is recommended as a positive control antibody for Western Blot analysis of enhanced human CRTAP expression in CRTAP transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

DATA



CRTAP (E-1): sc-393136. Western blot analysis of CRTAP expression in non-transfected: sc-117752 (A) and human CRTAP transfected: sc-170519 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.