



Legumain (h4): 293T Lysate: sc-171012

BACKGROUND

Legumain, also known as LGMN, AEP (asparaginyl endopeptidase) or PRSC1, is a 433 amino acid protein that localizes to the lysosome and belongs to the peptidase C13 family. Expressed ubiquitously with particularly high expression in placenta, heart and kidney, Legumain functions as a cysteine protease that specifically catalyzes the hydrolysis of asparaginyl and aspartyl bonds. Additionally, Legumain is thought to be involved in the processing of bacterial proteins for MHC class II antigen presentation in the lysosomal/endosomal system. Legumain exists as both a precursor and a fully mature, active enzyme that is produced in dendritic cells. Overexpression of Legumain may be associated with the formation of solid tumors, suggesting a role for Legumain in carcinogenesis. Multiple isoforms of Legumain exist due to alternative splicing events.

REFERENCES

1. Tanaka, T., Inazawa, J. and Nakamura, Y. 1996. Molecular cloning of a human cDNA encoding putative cysteine protease (PRSC1) and its chromosome assignment to 14q32.1. *Cytogenet. Cell Genet.* 74: 120-123.
2. Chen, J.M., Dando, P.M., Rawlings, N.D., Brown, M.A., Young, N.E., Stevens, R.A., Hewitt, E., Watts, C. and Barrett, A.J. 1997. Cloning, isolation, and characterization of mammalian Legumain, an asparaginyl endopeptidase. *J. Biol. Chem.* 272: 8090-8098.
3. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 602620. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
4. Li, D.N., Matthews, S.P., Antoniou, A.N., Mazzeo, D. and Watts, C. 2003. Multistep autoactivation of asparaginyl endopeptidase *in vitro* and *in vivo*. *J. Biol. Chem.* 278: 38980-38990.
5. Burster, T., Beck, A., Tolosa, E., Marin-Esteban, V., Röttschke, O., Falk, K., Lautwein, A., Reich, M., Brandenburg, J., Schwarz, G., Wiendl, H., Melms, A., Lehmann, R., Stevanovic, S., Kalbacher, H. and Driessen, C. 2004. Cathepsin G, and not the asparagine-specific endoprotease, controls the processing of myelin basic protein in lysosomes from human B lymphocytes. *J. Immunol.* 172: 5495-5503.
6. Murthy, R.V., Arberman, G., Gao, J., Roodman, G.D. and Sun, X.F. 2005. Legumain expression in relation to clinicopathologic and biological variables in colorectal cancer. *Clin. Cancer Res.* 11: 2293-2299.
7. Liu, Z., Jang, S.W., Liu, X., Cheng, D., Peng, J., Yepes, M., Li, X.J., Matthews, S., Watts, C., Asano, M., Hara-Nishimura, I., Luo, H.R. and Ye, K. 2008. Neuroprotective actions of PIKE-L by inhibition of SET proteolytic degradation by asparagine endopeptidase. *Mol. Cell* 29: 665-678.

CHROMOSOMAL LOCATION

Genetic locus: LGMN (human) mapping to 14q32.12.

PRODUCT

Legumain (h4): 293T Lysate represents a lysate of human Legumain transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

Legumain (h4): 293T Lysate is suitable as a Western Blotting positive control for human reactive Legumain antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.