

cyclin D1 (h4): 293T Lysate: sc-171254

BACKGROUND

The proliferation of eukaryotic cells is controlled at specific points in the cell cycle, particularly at the G₁ to S and the G₂ to M transitions. It is well established that the Cdc2 p34-cyclin B protein kinase plays a critical role in the G₂ to M transition, while cyclin A associates with Cdk2 p33 and functions in S phase. Considerable effort directed towards the identification of G₁ cyclins has led to the isolation of cyclin D, cyclin C and cyclin E. Of these, cyclin D corresponds to a putative human oncogene, designated PRAD1, which maps at the site of the Bcl-1 rearrangement in certain lymphomas and leukemias. Two additional human type D cyclins, as well as their mouse homologs, have been identified. Evidence has established that members of the cyclin D family function to regulate phosphorylation of the retinoblastoma gene product, thereby activating E2F transcription factors.

REFERENCES

1. Draetta, G. 1990. Cell cycle control in eukaryotes: molecular mechanisms of Cdc2 activation. *Trends Biol. Sci.* 15: 378-383.
2. Xiong, Y., et al. 1991. Human D-type cyclin. *Cell* 65: 691-699.
3. Inaba, T., et al. 1992. Genomic organization, chromosomal localization, and independent expression of human cyclin D genes. *Genomics* 13: 565-574.
4. Xiong, Y., et al. 1992. Molecular cloning and chromosomal mapping of CCND genes encoding human D-type cyclins. *Genomics* 13: 575-584.
5. Kiyokawa, H., et al. 1992. Cloning of a D-type cyclin from murine erythroleukemia cells. *Proc. Natl. Acad. Sci. USA* 89: 2444-2447.
6. Won, K., et al. 1992. Growth-regulated expression of D-type cyclin genes in human diploid fibroblasts. *Proc. Natl. Acad. Sci. USA* 89: 9910-9914.
7. Motokura, T., et al. 1992. Cloning and characterization of human cyclin D3, a cDNA closely related in sequence to the PRAD1/cyclin D1 proto-oncogene. *J. Biol. Chem.* 267: 20412-20415.
8. Ewen, M.E., et al. 1993. Functional interactions of the retinoblastoma protein with mammalian D-type cyclins. *Cell* 73: 487-497.
9. Dowdy, S.F., et al. 1993. Physical interaction of the retinoblastoma protein with human D cyclins. *Cell* 73: 499-511.

CHROMOSOMAL LOCATION

Genetic locus: CCND1 (human) mapping to 11q13.3.

PRODUCT

cyclin D1 (h4): 293T Lysate represents a lysate of human cyclin D1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

APPLICATIONS

cyclin D1 (h4): 293T Lysate is suitable as a Western Blotting positive control for human reactive cyclin D1 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.