



ZBED1 (h3): 293T Lysate: sc-172080

BACKGROUND

ZBED1 (zinc finger BED domain-containing protein 1), also known as ALTE (Ac-like transposable element), DREF or TRAMP, is a 694 amino acid protein that localizes specifically to granular structures within the nucleus. Expressed ubiquitously at low levels and present at higher levels in heart, placenta, spleen and skeletal muscle, ZBED1 is thought to function as a transcription factor that regulates a number of ribosomal protein (RP) encoding genes, thereby playing a role in the cell cycle and in cell proliferation events. ZBED1 contains one BED-type zinc finger and binds specifically to 5'-TGTCG[CT]GA [CT] A-3' DNA regions found in RP promoters. Additionally, ZBED1 binds strongly to the promoter region of Histone H1 (a protein required for the condensation of nucleosomes into higher order structures), subsequently activating H1 transcription.

REFERENCES

1. Oosumi, T., Belknap, W.R. and Garlick, B. 1995. Mariner transposons in humans. *Nature* 378: 672.
2. Nagase, T., Ishikawa, K., Suyama, M., Kikuno, R., Miyajima, N., Tanaka, A., Kotani, H., Nomura, N. and Ohara, O. 1998. Prediction of the coding sequences of unidentified human genes. XI. The complete sequences of 100 new cDNA clones from brain which code for large proteins *in vitro*. *DNA Res.* 5: 277-286.
3. Esposito, T., Gianfrancesco, F., Ciccodicola, A., Montanini, L., Mumm, S., D'Urso, M. and Forabosco, A. 1999. A novel pseudoautosomal human gene encodes a putative protein similar to Ac-like transposases. *Hum. Mol. Genet.* 8: 61-67.
4. Ohshima, N., Takahashi, M. and Hirose, F. 2003. Identification of a human homologue of the DREF transcription factor with a potential role in regulation of the Histone H1 gene. *J. Biol. Chem.* 278: 22928-22938.
5. Online Mendelian Inheritance in Man, OMIM™. 2003. Johns Hopkins University, Baltimore, MD. MIM Number: 300178. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
6. Yamashita, D., Komori, H., Higuchi, Y., Yamaguchi, T., Osumi, T. and Hirose, F. 2007. Human DNA replication-related element binding factor (hDREF) self-association via hATC domain is necessary for its nuclear accumulation and DNA binding. *J. Biol. Chem.* 282: 7563-7575.
7. Yamashita, D., Sano, Y., Adachi, Y., Okamoto, Y., Osada, H., Takahashi, T., Yamaguchi, T., Osumi, T. and Hirose, F. 2007. hDREF regulates cell proliferation and expression of ribosomal protein genes. *Mol. Cell. Biol.* 27: 2003-2013.

CHROMOSOMAL LOCATION

Genetic locus: ZBED1 (human) mapping to Xp22.33/Yp11.31.

PRODUCT

ZBED1 (h3): 293T Lysate represents a lysate of human ZBED1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

APPLICATIONS

ZBED1 (h3): 293T Lysate is suitable as a Western Blotting positive control for human reactive ZBED1 antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.