

PP2C α (h2): 293T Lysate: sc-172768

BACKGROUND

The phosphorylation and dephosphorylation of proteins on serine and threonine residues is an essential means of regulating a broad range of cellular functions in eukaryotes, including division, homeostasis and apoptosis. A group of proteins that are intimately involved in this process are the serine/threonine protein phosphatases. The protein phosphatase 2C α (PP2C α) has broad specificity. It dephosphorylates and negatively regulates the activities of MAP kinases and MAP kinase-kinases, and also inhibits the activation of p38 and JNK kinase cascades induced by environmental stresses. PP2C α also induces the expression of endogenous p53 and the p53-responsive gene p21, leading to cell cycle arrest and apoptosis. The PP2C α protein, which contains an active site containing a dinuclear metal ion center, shows highest expression in epithelial cells as well as in the digestive tract, lung, kidney, breast, prostate, endocrine glands and brain.

REFERENCES

1. Ueki, K., Muramatsu, T. and Kincaid, R.L. 1992. Structure and expression of two isoforms of the murine calmodulin-dependent protein phosphatase regulatory subunit (calcineurin B). *Biochem. Biophys. Res. Commun.* 187: 537-543.
2. Cohen, P.T. 1993. Important roles for novel protein phosphatases dephosphorylating serine and threonine residues. *Biochem. Soc. Trans.* 21: 884-888.
3. Yokoyama, N., Kobayashi, T., Tamura, S. and Sugiya, H. 1996. Purification and characterization of protein phosphatase 2C in rat parotid acinar cells: two forms of Mg²⁺-activated histone phosphatase and phosphorylation by cAMP-dependent protein kinase. *Arch. Biochem. Biophys.* 331: 1-8.
4. Takekawa, M., Maeda, T. and Saito, H. 1998. Protein phosphatase 2C α inhibits the human stress-responsive p38 and JNK MAPK pathways. *EMBO J.* 17: 4744-4452.
5. Lifschitz-Mercer, B., Sheinin, Y., Ben-Meir, D., Bramante-Schreiber, L., Leider-Trejo, L., Karby, S., Smorodinsky, N.I. and Lavi, S. 2001. Protein phosphatase 2C α expression in normal human tissues: an immunohistochemical study. *Histochem. Cell Biol.* 116: 31-39.
6. Jackson, M.D., Fjeld, C.C. and Denu, J.M. 2003. Probing the function of conserved residues in the serine/threonine phosphatase PP2C α . *Biochemistry* 42: 8513-8521.
7. Ofek, P., Ben-Meir, D., Kariv-Inbal, Z., Oren, M. and Lavi, S. 2003. Cell cycle regulation and p53 activation by protein phosphatase 2C α . *J. Biol. Chem.* 278: 14299-14305.

CHROMOSOMAL LOCATION

Genetic locus: PPM1A (human) mapping to 14q23.1.

PRODUCT

PP2C α (h2): 293T Lysate represents a lysate of human PP2C α transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

RESEARCH USE

For research use only, not for use in diagnostic procedures.


APPLICATIONS

PP2C α (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive PP2C α antibodies. Recommended use: 10-20 μ l per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

PP2C α (p6c7): sc-56956 is recommended as a positive control antibody for Western Blot analysis of enhanced human PP2C α expression in PP2C α transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

DATA



PP2C α (p6c7): sc-56956. Western blot analysis of PP2C α expression in non-transfected: sc-117752 (A) and human PP2C α transfected: sc-172768 (B) 293T whole cell lysates.



PP2C α (6D708): sc-71922. Western blot analysis of PP2C α expression in non-transfected: sc-117752 (A) and human PP2C α transfected: sc-172768 (B) 293T whole cell lysates.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.