

Cdc5L (h2): 293T Lysate: sc-173633

BACKGROUND

Cdc5L (cell division cycle 5-like protein, *Pombe* Cdc5-related protein) is a DNA-binding protein encoded by the human gene CDC5L. Cdc5L contains two HTH Myb-type DNA-binding domains and may shuttle between cytoplasm and nucleus. It is involved in cell cycle control and may act as a transcription activator. Cdc5L is a spliceosomal protein that is highly conserved across species. It belongs to the group of proteins that comprise the core of spliceosomal complexes and are essential for pre-mRNA splicing. Cdc5L is involved in the second catalytic step of pre-mRNA splicing, which involves cleavage at the 3' splice site and the ligation of the exons. This process releases the intact intron lariat. A chromosomal aberration involving Cdc5L is found in multicystic renal dysplasia. This aberration is caused by a translocation (t(6;19;p21;q13.1)) with USF-2.

REFERENCES

1. Hirayama, T. and Shinozaki, K. 1996. A Cdc5⁺ homolog of a higher plant, *Arabidopsis thaliana*. Proc. Natl. Acad. Sci. USA 93: 13371-13376.
2. Bernstein, H.S. and Coughlin, S.R. 1997. *Pombe* Cdc5-related protein. A putative human transcription factor implicated in mitogen-activated signaling. J. Biol. Chem. 272: 5833-5837.
3. Groenen, P.M., Vanderlinden, G., Devriendt, K., Fryns, J.P. and Van de Ven, W.J. 1998. Rearrangement of the human CDC5L gene by a t(6;19)(p21;q13.1) in a patient with multicystic renal dysplasia. Genomics 49: 218-229.
4. Ajuh, P., Kuster, B., Panov, K., Zomerdijk, J.C., Mann, M. and Lamond, A.I. 2000. Functional analysis of the human Cdc5L complex and identification of its components by mass spectrometry. EMBO J. 19: 6569-6581.
5. Ajuh, P., Sleeman, J., Chusainow, J. and Lamond, A.I. 2001. A direct interaction between the carboxyl-terminal region of Cdc5L and the WD40 domain of PLRG1 is essential for pre-mRNA splicing. J. Biol. Chem. 276: 42370-42381.
6. Jurica, M.S., Licklider, L.J., Gygi, S.R., Grigorieff, N. and Moore, M.J. 2002. Purification and characterization of native spliceosomes suitable for three-dimensional structural analysis. RNA 8: 426-439.
7. Leonard, D., Ajuh, P., Lamond, A.I. and Legerski, R.J. 2003. hLodestar/HuF2 interacts with Cdc5L and is involved in pre-mRNA splicing. Biochem. Biophys. Res. Commun. 308: 793-801.

CHROMOSOMAL LOCATION

Genetic locus: CDC5L (human) mapping to 6p21.1.

PRODUCT

Cdc5L (h2): 293T Lysate represents a lysate of human Cdc5L transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.


APPLICATIONS

Cdc5L (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive Cdc5L antibodies. Recommended use: 10-20 µl per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

Cdc5L (2136C1a): sc-81220 is recommended as a positive control antibody for Western Blot analysis of enhanced human Cdc5L expression in Cdc5L transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

DATA



Cdc5L (2136C1a): sc-81220. Western blot analysis of Cdc5L expression in non-transfected: sc-117752 (A) and human truncated Cdc5L transfected: sc-173633 (B) 293T whole cell lysates.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.