



AKAP 7 (h): 293T Lysate: sc-174003

BACKGROUND

The type II cAMP-protein kinase (PKA) is a multifunctional kinase with a broad range of substrates. Specificity of PKA signaling is thought to be mediated by the compartmentalization of the kinase to specific sites within the cell. To maintain this specific localization, the regulatory (R) subunits (RI and RII) of PKA interact with specific R-anchoring proteins designated AKAPs (A-kinase anchoring proteins). AKAP 7 (A-kinase anchor protein 7), also known as AKAP18, is a 104 amino acid protein that belongs to the AKAP family. AKAP 7 is expressed in brain, heart, lung, pancreas and skeletal muscle. AKAP 7 binds PKA to the plasma membrane, and permits functional coupling to the L-type calcium channel. Four isoforms exist due to alternative splicing events. It has been suggested that the γ isoform binds RI and may be responsible for positioning PKA via RI and/or RII to regulate PKA-mediated gene transcription in both somatic cells and oocytes.

REFERENCES

1. Trotter, K.W., et al. 1999. Alternative splicing regulates the subcellular localization of A-kinase anchoring protein 18 isoforms. *J. Cell Biol.* 147: 1481-1492.
2. Scott, J.D., et al. 2000. Coordination of cAMP signaling events through PKA anchoring. *Adv. Pharmacol.* 47: 175-207.
3. Edwards, A.S., et al. 2000. A-kinase anchoring proteins: protein kinase A and beyond. *Curr. Opin. Cell Biol.* 12: 217-221.
4. Klusmann, E., et al. 2001. Role and identification of protein kinase A anchoring proteins in vasopressin-mediated aquaporin-2 translocation. *Kidney Int.* 60: 446-449.
5. Brown, R.L., et al. 2003. AKAP7 γ is a nuclear RI-binding AKAP. *Biochem. Biophys. Res. Commun.* 306: 394-401.
6. Henn, V., et al. 2004. Identification of a novel A-kinase anchoring protein 18 isoform and evidence for its role in the vasopressin-induced aquaporin-2 shuttle in renal principal cells. *J. Biol. Chem.* 279: 26654-26665.
7. Hundsruker, C., et al. 2006. High-affinity AKAP7 δ -protein kinase A interaction yields novel protein kinase A-anchoring disruptor peptides. *Biochem. J.* 396: 297-306.
8. McSorley, T., et al. 2006. Spatial organisation of AKAP18 and PDE4 isoforms in renal collecting duct principal cells. *Eur. J. Cell Biol.* 85: 673-678.
9. Gold, M.G., et al. 2008. AKAP18 contains a phosphoesterase domain that binds AMP. *J. Mol. Biol.* 375: 1329-1343.

CHROMOSOMAL LOCATION

Genetic locus: AKAP7 (human) mapping to 6q23.2.

PRODUCT

AKAP 7 (h): 293T Lysate represents a lysate of human AKAP 7 transfected 293T cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

APPLICATIONS

AKAP 7 (h): 293T Lysate is suitable as a Western Blotting positive control for human reactive AKAP 7 antibodies. Recommended use: 10-20 μ l per lane.

Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.