

# 4E-BP1 (h2): 293T Lysate: sc-174963

## BACKGROUND

The translation of proteins from eukaryotic mRNA is initiated by the multi-subunit complex eIF-4F, which associates with the mRNA 5' cap structure. eIF-4E, a component of eIF-4F, is responsible for binding to the 5' cap structure and for the assembly of the eIF-4F complex. The regulatory protein 4E-BP1, also referred to as PHAS-I, inhibits eIF-4E function. Phosphorylation of 4E-BP1 by S6 kinase p70, MAP kinases or PKCs causes the disassociation of 4E-BP1 from eIF-4E, promoting translation. A protein that is functionally related to 4E-BP1, designated 4E-BP2, also associates with eIF-4E.

## REFERENCES

- Lin, T.A., Kong, X., Haystead, T.A., Pause, A., Belsham, G., Sonenberg, N. and Lawrence, J.C., Jr. 1994. PHAS-I as a link between mitogen-activated protein kinase and translation initiation. *Science* 266: 653-656.
- Rau, M., Ohlmann, T., Morley, S.J. and Pain, V.M. 1996. A reevaluation of the cap-binding protein, eIF-4E, as a rate-limiting factor for initiation of translation in reticulocyte lysate. *J. Biol. Chem.* 271: 8983-8990.
- Whalen, S.G., Gingras, A.C., Amankwa, L., Mader, S., Branton, P.E., Aebersold, R. and Sonenberg, N. 1996. Phosphorylation of eIF-4E on serine 209 by protein kinase C is inhibited by the translational repressors, 4E-binding proteins. *J. Biol. Chem.* 271: 11831-11837.
- Diggle, T.A., Moule, S.K., Avison, M.B., Flynn, A., Foulstone, E.J., Proud, C.G. and Denton, R.M. 1996. Both Rapamycin-sensitive and -insensitive pathways are involved in the phosphorylation of the initiation factor-4E-binding protein (4E-BP1) in response to Insulin in rat epididymal fat-cells. *Biochem. J.* 316: 447-453.
- Beretta, L., Gingras, A.C., Svitkin, Y.V., Hall, M.N. and Sonenberg, N. 1996. Rapamycin blocks the phosphorylation of 4E-BP1 and inhibits cap-dependent initiation of translation. *EMBO J.* 15: 658-664.
- Mendez, R., Myers, M.G., Jr., White, M.F. and Rhoads, R.E. 1996. Stimulation of protein synthesis, eukaryotic translation initiation factor 4E phosphorylation, and PHAS-I phosphorylation by Insulin requires Insulin receptor substrate 1 and phosphatidylinositol 3-kinase. *Mol. Cell. Biol.* 16: 2857-2864.
- von Manteuffel, S.R., Gingras, A.C., Ming, X.F., Sonenberg, N. and Thomas, G. 1996. 4E-BP1 phosphorylation is mediated by the FRAP-p70s6k pathway and is independent of mitogen-activated protein kinase. *Proc. Natl. Acad. Sci. USA* 93: 4076-4080.

## CHROMOSOMAL LOCATION

Genetic locus: EIF4EBP1 (human) mapping to 8p11.23.

## PRODUCT

4E-BP1 (h2): 293T Lysate represents a lysate of human 4E-BP1 transfected 293T cells and is provided as 100 µg protein in 200 µl SDS-PAGE buffer.

## STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

## APPLICATIONS

4E-BP1 (h2): 293T Lysate is suitable as a Western Blotting positive control for human reactive 4E-BP1 antibodies. Recommended use: 10-20 µl per lane.


Control 293T Lysate: sc-117752 is available as a Western Blotting negative control lysate derived from non-transfected 293T cells.

4E-BP1 (11G12C11): sc-81149 is recommended as a positive control antibody for Western Blot analysis of enhanced human 4E-BP1 expression in 4E-BP1 transfected 293T cells (starting dilution 1:100, dilution range 1:100-1:1,000).

## RECOMMENDED SUPPORT REAGENTS

To ensure optimal results, the following support reagents are recommended:  
1) Western Blotting: use m-IgGκ BP-HRP: sc-516102 or m-IgGκ BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048.

## DATA



4E-BP1 (11G12C11): sc-81149. Western blot analysis of 4E-BP1 expression in non-transfected: sc-117752 (A) and human 4E-BP1 transfected: sc-174963 (B) 293T whole cell lysates.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.