

TAAR4 (m): 293 Lysate: sc-179563

BACKGROUND

Trace amines are endogenous molecules structurally related to classical biogenic amines that are linked to psychiatric conditions. TAAR4 belongs to a family of G protein-coupled receptors, referred to as trace-amine-associated receptors (TAARs), which are activated by trace amines and are present in very low levels in mammalian tissue. TAARs contain several structural features that are similar to the rhodopsin β -adrenergic receptor superfamily, including the positions of the seven-transmembrane regions that provide common ligand-binding pockets as well as the short N- and C-terminal domains. TAAR proteins are potential targets for studying amine-containing drugs of abuse, such as amphetamines and MDMA, as well as neuropsychiatric disorders including schizophrenia, depression and attention deficit disorder.

REFERENCES

1. Online Mendelian Inheritance in Man, OMIM™. 2002. Johns Hopkins University, Baltimore, MD. MIM Number: 604849. World Wide Web URL: <http://www.ncbi.nlm.nih.gov/omim/>
2. Aridon, P., D'Andrea, G., Rigamonti, A., Leone, M., Casari, G. and Bussone, G. 2004. Elusive amines and cluster headache: mutational analysis of trace amine receptor cluster on chromosome 6q23. *Neurol. Sci.* 3: S279-S280.
3. Lindemann, L., Ebeling, M., Kratochwil, N.A., Bunzow, J.R., Grandy, D.K. and Hoener, M.C. 2005. Trace amine-associated receptors form structurally and functionally distinct subfamilies of novel G protein-coupled receptors. *Genomics* 85: 372-385.
4. Miller, G.M., Verrico, C.D., Jassen, A., Konar, M., Yang, H., Panas, H., Bahn, M., Johnson, R. and Madras, B.K. 2005. Primate trace amine receptor 1 modulation by the dopamine transporter. *J. Pharmacol. Exp. Ther.* 313: 983-994.
5. Liberles, S.D. and Buck, L.B. 2006. A second class of chemosensory receptors in the olfactory epithelium. *Nature* 442: 645-650.
6. Wainscott, D.B., Little, S.P., Yin, T., Tu, Y., Rocco, V.P., He, J.X. and Nelson, D.L. 2007. Pharmacologic characterization of the cloned human trace amine-associated receptor 1 (TAAR1) and evidence for species differences with the rat TAAR1. *J. Pharmacol. Exp. Ther.* 320: 475-485.

CHROMOSOMAL LOCATION

Genetic locus: Taar4 (mouse) mapping to 10 A4.

PRODUCT

TAAR4 (m): 293 Lysate represents a lysate of mouse TAAR4 transfected 293 cells and is provided as 100 μ g protein in 200 μ l SDS-PAGE buffer.

STORAGE

Store at -20° C. Repeated freezing and thawing should be minimized. Sample vial should be boiled once prior to use. Non-hazardous. No MSDS required.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

APPLICATIONS

TAAR4 (m): 293 Lysate is suitable as a Western Blotting positive control for mouse reactive TAAR4 antibodies. Recommended use: 10-20 μ l per lane.

Control 293 Lysate: sc-110760 is available as a Western Blotting negative control lysate derived from non-transfected 293 cells.

RESEARCH USE

For research use only, not for use in diagnostic procedures.