

normal mouse IgG: sc-2025

BACKGROUND

Santa Cruz Biotechnology offers a wide variety of control immunoglobulin and control sera for a large selection of species, including mouse, rabbit, goat, chicken, rat, hamster, canine, guinea pig and ovine. Control immunoglobulin and immunoglobulin conjugates are useful negative controls. Normal sera is offered to be used as blocking reagents. Santa Cruz Biotechnology offers affinity purified normal immunoglobulins and immunoglobulin conjugates for use as negative controls in applications including flow cytometry, immunohistochemistry and immunofluorescence. Isotype specific control immunoglobulins include classes such as mouse IgG₁, IgG_{2a}, IgG_{2b}, IgG₃, IgM and IgA, rat IgG₁, IgG_{2a}, IgG_{2b} and IgM, Armenian hamster IgG, and both goat and rabbit IgG. Most control immunoglobulins are available as unconjugated controls or as FITC (fluorescein isothiocyanate), PE (phycoerythrin), PE-Cy5 (phycoerythrin-Cy5), PE-Cy7 (phycoerythrin-Cy7), APC (allophycocyanin) and APC-Cy7 (allophycocyanin-Cy7) conjugates. Additional conjugates include Alexa Fluor[®] 488, Alexa Fluor[®] 647, Alexa Fluor[®] 405, PerCP (peridinin chlorophyll protein complex) and PerCP-Cy5.5 (peridinin chlorophyll protein complex-Cy 5.5).

SOURCE

normal mouse IgG is an affinity purified, unconjugated conjugated isotype control immunoglobulin from mouse.

PRODUCT

Each vial contains 200 µg mouse IgG in PBS containing 0.1% gelatin and 0.1% sodium azide.

APPLICATIONS

normal mouse IgG is recommended for use as an isotype control immunoglobulin in place of a target specific primary antibody of the same isotype (mouse IgG) by Western Blotting, immunoprecipitation and immunohistochemical staining. To be used at an assay dependent dilution.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

SELECT PRODUCT CITATIONS

1. Tan, M., et al. 2002. Phosphorylation on tyrosine-15 of p34(Cdc2) by ErbB2 inhibits p34(Cdc2) activation and is involved in resistance to taxol-induced apoptosis. *Mol. Cell* 9: 993-1004.
2. Acevedo-Duncan, M., et al. 2002. Human glioma PKC- ι and PKC- β II phosphorylate cyclin-dependent kinase activating kinase during the cell cycle. *Cell Prolif.* 35: 23-36.
3. Xie, R., et al. 2011. Microtubule-associated protein 1S (MAP1S) bridges autophagic components with microtubules and mitochondria to affect autophagosomal biogenesis and degradation. *J. Biol. Chem.* 286: 10367-10377.
4. Bernatik, O., et al. 2011. Sequential activation and inactivation of dishevelled in the Wnt/ β -catenin pathway by casein kinases. *J. Biol. Chem.* 286: 10396-10410.
5. Lenasi, T., et al. 2011. Cap-binding protein complex links pre-mRNA capping to transcription elongation and alternative splicing through positive transcription elongation factor β (P-TEF β). *J. Biol. Chem.* 286: 22758-22768.
6. Gizard, F., et al. 2011. Telomerase activation in atherosclerosis and induction of telomerase reverse transcriptase expression by inflammatory stimuli in macrophages. *Arterioscler. Thromb. Vasc. Biol.* 31: 245-252.
7. Euskirchen, G.M., et al. 2011. Diverse roles and interactions of the SWI/SNF chromatin remodeling complex revealed using global approaches. *PLoS Genet.* 7: e1002008.
8. Li, B., et al. 2011. Increased hedgehog signaling in postnatal kidney results in aberrant activation of nephron developmental programs. *Hum. Mol. Genet.* 20: 4155-4166.
9. Liu, J., et al. 2011. Hyperglycemia-induced cerebral hematoma expansion is mediated by plasma kallikrein. *Nat. Med.* 17: 206-210.
10. Calaf, G.M., et al. 2011. Protective role of curcumin in oxidative stress of breast cells. *Oncol. Rep.* 26: 1029-1035.
11. Song, Y., et al. 2011. Ligand-dependent corepressor acts as a novel corepressor of thyroid hormone receptor and represses hepatic lipogenesis in mice. *J. Hepatol.* 56: 248-254.
12. Taura, M., et al. 2011. MEF/ELF4 transactivation by E2F1 is inhibited by p53. *Nucleic Acids Res.* 39: 76-88.
13. Chen, S.T., et al. 2012. Recombinant MPT83 derived from *Mycobacterium tuberculosis* induces cytokine production and upregulates the function of mouse macrophages through TLR2. *J. Immunol.* 188: 668-677.
14. Reid, H.M., et al. 2012. Interaction of the human prostacyclin receptor and the NHERF4 family member intestinal and kidney enriched PDZ protein (IKEPP). *Biochim. Biophys. Acta* 1823: 1998-2012.
15. Zanin-Zhorov, A., et al. 2012. Scaffold protein Disc large homolog 1 is required for T-cell receptor-induced activation of regulatory T-cell function. *Proc. Natl. Acad. Sci. USA* 109: 1625-1630.
16. Kim, Y.C., et al. 2012. Angiotensin II regulates activation of Bim via Rb/E2F1 during apoptosis: involvement of interaction between AMPK β 1/2 and Cdk4. *Am. J. Physiol. Lung Cell Mol. Physiol.* 303: L228-L238.