

normal rabbit IgG: sc-2027

BACKGROUND

Santa Cruz Biotechnology offers a wide variety of control immunoglobulin and control sera for a large selection of species, including mouse, rabbit, goat, chicken, rat, hamster, canine, guinea pig and sheep. Control immunoglobulin and immunoglobulin conjugates are useful negative controls. Normal sera is offered to be used as blocking reagents. Santa Cruz Biotechnology offers affinity purified normal immunoglobulins and immunoglobulin conjugates for use as negative controls in applications including flow cytometry, immunohistochemistry and immunofluorescence. Isotype specific control immunoglobulins include classes such as mouse IgG₁, IgG_{2a}, IgG_{2b}, IgG₃, IgM and IgA, rat IgG₁, IgG_{2a}, IgG_{2b} and IgM, Armenian hamster IgG, and both goat and rabbit IgG. Most control immunoglobulins are available as unconjugated controls or as FITC (fluorescein isothiocyanate), PE (phycoerythrin), PE-Cy5 (phycoerythrin-Cy5), PE-Cy7 (phycoerythrin-Cy7), APC (allophycocyanin) and APC-Cy7 (allophycocyanin-Cy7) conjugates. Additional conjugates include Alexa Fluor[®] 488, Alexa Fluor[®] 647, Alexa Fluor[®] 405, PerCP (peridinin chlorophyll protein complex) and PerCP-Cy5.5 (peridinin chlorophyll protein complex-Cy 5.5).

SOURCE

normal rabbit IgG is an unconjugated, affinity purified isotype control immunoglobulin from rabbit.

PRODUCT

Each vial contains 200 µg rabbit IgG in 0.5 ml PBS with 0.1% sodium azide and 0.1% gelatin.

APPLICATIONS

normal rabbit IgG is recommended for use as a isotype control immunoglobulin in place of a target specific primary antibody of the same isotype (rabbit IgG) by Western blotting, immunofluorescence staining, immunohistochemical staining and flow cytometry. To be used at an assay dependent dilution.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.

SELECT PRODUCT CITATIONS

- Mal, A., et al. 2001. A role for histone deacetylase HDAC1 in modulating the transcriptional activity of MyoD: inhibition of the myogenic program. *EMBO J.* 20: 1739-1753.
- Naren, A.P., et al. 2003. A macromolecular complex of β2 adrenergic receptor, CFTR, and ezrin/radixin/moesin-binding phosphoprotein 50 is regulated by PKA. *Proc. Natl. Acad. Sci. USA* 100: 342-346.
- Sepulveda, M.A., et al. 2004. NFκB and Oct-2 synergize to activate the human 3' Igh hs4 enhancer in B cells. *J. Immunol.* 172: 1054-1064.
- Herrera, V.L., et al. 2005. Embryonic lethality in *Dear* gene-deficient mice: new player in angiogenesis. *Am. J. Physiol.* 23: 257-268.
- Deschênes, J. and Bourdeau, V. 2007. Regulation of GREB1 transcription by estrogen receptor α through a multipartite enhancer spread over 20 kb of upstream flanking sequences. *J. Biol. Chem.* 282: 17335-17339.
- Liu, H., et al. 2008. The FGL2-FcγRIIB pathway: a novel mechanism leading to immunosuppression. *Eur. J. Immunol.* 38: 3114-3126.
- Rouget, C., et al. 2010. Maternal mRNA deadenylation and decay by the πRNA pathway in the early *Drosophila* embryo. *Nature* 467: 1128-1132.
- Udo, K., et al. 2010. Adipose tissue explants and MDCK cells reciprocally regulate their morphogenesis in coculture. *Kidney Int.* 78: 60-68.
- Huang, L., et al. 2011. Prevention of transcriptional silencing by a replicator-binding complex consisting of SWI/SNF, MeCP1, and hnRNP C1/C2. *Mol. Cell. Biol.* 31: 3472-3484.
- Koch, D., et al. 2011. Proper synaptic vesicle formation and neuronal network activity critically rely on syndapin I. *EMBO J.* 30: 4955-4969.
- Kim, J.H., et al. 2011. Brain-derived neurotrophic factor uses CREB and Egr3 to regulate NMDA receptor levels in cortical neurons. *J. Neurochem.* 120: 210-219.
- Alimirah, F., et al. 2011. Functional significance of vitamin D receptor FokI polymorphism in human breast cancer cells. *PLoS ONE* 6: e16024.
- Strakovsky, R.S., et al. 2011. Gestational high fat diet programs hepatic phosphoenolpyruvate carboxykinase gene expression and histone modification in neonatal offspring rats. *J. Physiol.* 589: 2707-2717.
- Alberi, L., et al. 2011. Activity-induced Notch signaling in neurons requires Arc/Arg3.1 and is essential for synaptic plasticity in hippocampal networks. *Neuron* 69: 437-444.
- Keating, G., et al. 2012. Regulation of the human prostacyclin receptor gene in megakaryocytes: Major roles for C/EBPδ and PU.1. *Biochim Biophys Acta.* 1819: 428-445.
- Song, Y., et al. 2012. Ligand-dependent corepressor acts as a novel corepressor of thyroid hormone receptor and represses hepatic lipogenesis in mice. *J. Hepatol.* 56: 248-254.
- Zhang, X., et al. 2012. Hepatic cellular senescence pathway genes are induced through histone modifications in a diet-induced obese rat model. *Am. J. Physiol. Gastrointest. Liver Physiol.* 302: G558-G564.