

CNKSR3 (E-20): sc-242478

BACKGROUND

CNKSR3 (connector enhancer of kinase suppressor of ras 3), also known as CNK homolog protein 3, CNKSR family member 3 or maguin-like protein, is a 555 amino acid protein that belongs to the CNKSR family. Suggested to play a role in transepithelial sodium transport, CNKSR3 mediates aldosterone-induced and ENaC-mediated sodium transport, and is upregulated by aldosterone. CNKSR3 contains one SAM (sterile α motif) domain, a CRIC domain, DUF1170 domain and a single PDZ (DHR) domain. The gene encoding CNKSR3 maps to human chromosome 6, which contains 170 million base pairs and comprises nearly 6% of the human genome. Deletion of a portion of the q arm of chromosome 6 is associated with early onset intestinal cancer, suggesting the presence of a cancer susceptibility locus. Additionally, Porphyria cutanea tarda, Parkinson's disease, Stickler syndrome and a susceptibility to bipolar disorder are all associated with genes that map to chromosome 6.

REFERENCES

1. Brunner, H.G., van Beersum, S.E., Warman, M.L., Olsen, B.R., Ropers, H.H. and Mariman, E.C. 1994. A Stickler syndrome gene is linked to chromosome 6 near the COL11A2 gene. *Hum. Mol. Genet.* 3: 1561-1564.
2. Cesari, R., Martin, E.S., Calin, G.A., Pentimalli, F., Bichi, R., McAdams, H., Trapasso, F., Drusco, A., Shimizu, M., Masciullo, V., D'Andrilli, G., Scambia, G., Picchio, M.C., Alder, H., Godwin, A.K. and Croce, C.M. 2003. Parkin, a gene implicated in autosomal recessive juvenile parkinsonism, is a candidate tumor suppressor gene on chromosome 6q25-q27. *Proc. Natl. Acad. Sci. USA* 100: 5956-5961.
3. Wright, G.J., Leslie, J.D., Ariza-McNaughton, L. and Lewis, J. 2004. Delta proteins and MAGI proteins: an interaction of Notch ligands with intracellular scaffolding molecules and its significance for zebrafish development. *Development* 131: 5659-5669.
4. Bläker, H., Mechttersheimer, G., Sutter, C., Hertkorn, C., Kern, M.A., Rieker, R.J., Penzel, R., Schirmacher, P. and Kloor, M. 2008. Recurrent deletions at 6q in early age of onset non-HNPCC- and non-FAP-associated intestinal carcinomas. Evidence for a novel cancer susceptibility locus at 6q14-q22. *Genes Chromosomes Cancer* 47: 159-164.
5. Ziera, T., Irlbacher, H., Fromm, A., Latouche, C., Krug, S.M., Fromm, M., Jaisser, F. and Borden, S.A. 2009. Cnksr3 is a direct mineralocorticoid receptor target gene and plays a key role in the regulation of the epithelial sodium channel. *FASEB J.* 23: 3936-3946.
6. Jalil, S., Grady, J.J., Lee, C. and Anderson, K.E. 2010. Associations among behavior-related susceptibility factors in porphyria cutanea tarda. *Clin. Gastroenterol. Hepatol.* 8: 297-302.

CHROMOSOMAL LOCATION

Genetic locus: CNKSR3 (human) mapping to 6q25.2; Cnksr3 (mouse) mapping to 10 A1.

SOURCE

CNKSR3 (E-20) is an affinity purified goat polyclonal antibody raised against a peptide mapping within an internal region of CNKSR3 of human origin.

PRODUCT

Each vial contains 200 μ g IgG in 1.0 ml of PBS with < 0.1% sodium azide and 0.1% gelatin.

Blocking peptide available for competition studies, sc-242478 P, (100 μ g peptide in 0.5 ml PBS containing < 0.1% sodium azide and 0.2% BSA).

APPLICATIONS

CNKSR3 (E-20) is recommended for detection of CNKSR3 of mouse, rat and human origin by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000), immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500) and solid phase ELISA (starting dilution 1:30, dilution range 1:30-1:3000).

CNKSR3 (E-20) is also recommended for detection of CNKSR3 in additional species, including equine, canine, bovine, porcine and avian.

Suitable for use as control antibody for CNKSR3 siRNA (h): sc-95360, CNKSR3 siRNA (m): sc-142435, CNKSR3 shRNA Plasmid (h): sc-95360-SH, CNKSR3 shRNA Plasmid (m): sc-142435-SH, CNKSR3 shRNA (h) Lentiviral Particles: sc-95360-V and CNKSR3 shRNA (m) Lentiviral Particles: sc-142435-V.

Molecular Weight of CNKSR3: 62 kDa.

RECOMMENDED SECONDARY REAGENTS

To ensure optimal results, the following support (secondary) reagents are recommended: 1) Western Blotting: use donkey anti-goat IgG-HRP: sc-2020 (dilution range: 1:2000-1:100,000) or Cruz Marker™ compatible donkey anti-goat IgG-HRP: sc-2033 (dilution range: 1:2000-1:5000), Cruz Marker™ Molecular Weight Standards: sc-2035, TBS Blotto A Blocking Reagent: sc-2333 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use donkey anti-goat IgG-FITC: sc-2024 (dilution range: 1:100-1:400) or donkey anti-goat IgG-TR: sc-2783 (dilution range: 1:100-1:400) with UltraCruz™ Mounting Medium: sc-24941.

STORAGE

Store at 4° C, ****DO NOT FREEZE****. Stable for one year from the date of shipment. Non-hazardous. No MSDS required.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com or our catalog for detailed protocols and support products.