



# CD2AP siRNA (r): sc-270133

## BACKGROUND

CD2-associated protein (CD2AP) is a cytosolic adaptor molecule that complexes with the intracellular portion of CD2 upon T cell activation. T cell activation induces cell adhesion through CD2-mediated binding to surface ligands on antigen-presenting cells, which enhances antigen-specific T cell activation, potentiates cell clustering and induces cytoskeletal polarization. CD2AP is expressed at highest levels in liver, thymus and spleen. CD2AP contains three SH3 domains that are essential for the interaction with CD2. Mutations in CD2AP that impair this interaction result in the disruption of cell clustering and polarization in activated T lymphocytes. Mice deficient in CD2AP develop a lethal congenital nephrotic syndrome, indicating that CD2AP is also involved in maintaining the integrity of the renal glomerulus.

## REFERENCES

1. Shaw, A.S., et al. 1997. Making the T cell receptor go the distance: a topological view of T cell activation. *Immunity* 6: 361-369.
2. Dustin, M.L., et al. 1998. A novel adaptor protein orchestrates receptor patterning and cytoskeletal polarity in T cell contacts. *Cell* 94: 667-677.
3. Nishizawa, K., et al. 1998. Identification of a proline-binding motif regulating CD2-triggered T lymphocyte activation. *Proc. Natl. Acad. Sci. USA* 95: 14897-14902.
4. Shih, N.Y., et al. 1999. Congenital nephrotic syndrome in mice lacking CD2-associated protein. *Science* 286: 312-315.
5. Guan, F., et al. 2006. Autocrine VEGF-A system in podocytes regulates podocin and its interaction with CD2AP. *Am. J. Physiol. Renal Physiol.* 291: F422-F428.
6. Fan, Q., et al. 2006. The relationship among nephrin, podocin, CD2AP, and  $\alpha$ -actinin might not be a true "interaction" in podocyte. *Kidney Int.* 69: 1207-1215.
7. Xia, W., et al. 2006. Interactions between transforming growth factor- $\beta$ 3/T $\beta$ R1, TAB1, and CD2AP disrupt blood-testis barrier and Sertoli-germ cell adhesion. *J. Biol. Chem.* 281: 16799-16813.
8. Konishi, H., et al. 2006. CFBP is a novel tyrosine-phosphorylated protein that might function as a regulator of CIN85/CD2AP. *J. Biol. Chem.* 281: 28919-28931.
9. Tossidou, I., et al. 2007. CD2AP/CIN85 balance determines receptor tyrosine kinase signaling response in podocytes. *J. Biol. Chem.* 282: 7457-7464.

## CHROMOSOMAL LOCATION

Genetic locus: Cd2ap (rat) mapping to 9q12.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## PRODUCT

CD2AP siRNA (r) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CD2AP shRNA Plasmid (r): sc-270133-SH and CD2AP shRNA (r) Lentiviral Particles: sc-270133-V as alternate gene silencing products.

For independent verification of CD2AP (r) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-270133A, sc-270133B and sc-270133C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CD2AP siRNA (r) is recommended for the inhibition of CD2AP expression in rat cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CD2AP gene expression knockdown using RT-PCR Primer: CD2AP (r)-PR: sc-270133-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.