

CyPB siRNA (h2): sc-270696

BACKGROUND

The immunosuppressant cyclosporin A (CsA) forms a trimolecular complex with cyclophilin and calcineurins to inhibit calcineurin phosphatase activity. Cyclophilins are conserved, ubiquitous and abundant cytosolic peptidyl-prolyl *cis-trans* isomerases that accelerate the isomerization of XaaPro peptide bonds and the refolding of proteins. Human cyclophilin A (CyPA), an intracellular protein of 165 amino acids, is the target of the Cyclosporin A and is encoded by a single unique gene conserved from yeast to humans. CyPA is known for its involvement in T cell differentiation and proliferation and is highly expressed in brain. CyPA is incorporated into the virion of the type 1 human immunodeficiency virus (HIV-1) via a direct interaction with the capsid domain of the viral Gag polyprotein and is crucial for efficient viral replication. Cyclophilin B (CyPB) is a member of the cyclophilin family with specific N- and C-terminal extensions. Unlike CyPA, CyPB has a signal sequence leading to its translocation in the endoplasmic reticulum. CyPB is secreted in biological fluids such as blood or milk and binds to a specific receptor present on the human lymphoblastic cell line Jurkat and on human peripheral blood lymphocytes.

REFERENCES

- Hasel, K.W., et al. 1991. An endoplasmic reticulum-specific cyclophilin. *Mol. Cell. Biol.* 11: 3484-3491.
- Arber, S., et al. 1992. s-cyclophilin is retained intracellularly via a unique COOH-terminal sequence and colocalizes with the calcium storage protein calreticulin. *J. Cell Biol.* 116: 113-125.
- Pflugl, G., et al. 1993. X-ray structure of a decameric cyclophilin-cyclosporin crystal complex. *Nature* 361: 91-94.
- Le Hir, M., et al. 1995. *In situ* detection of cyclosporin A: evidence for nuclear localization of cyclosporine and cyclophilins. *Lab. Invest.* 73: 727-733.
- Mariller, C., et al. 1996. Involvement of the N-terminal part of cyclophilin B in the interaction with specific Jurkat T-cell binding sites. *Biochem. J.* 317: 571-576.
- Mariller, C., et al. 1996. Evidence that human milk isolated cyclophilin B corresponds to a truncated form. *Biochim. Biophys. Acta* 1293: 31-38.
- Vajdos, F.F., et al. 1997. Crystal structure of cyclophilin A complexed with a binding site peptide from the HIV-1 capsid protein. *Protein Sci.* 6: 2297-2307.
- Endrich, M.M., et al. 1999. Maturation-induced conformational changes of HIV-1 capsid protein and identification of two highly affinity sites for cyclophilins in the C-terminal domain. *J. Biol. Chem.* 274: 5326-5332.
- Wang, P., et al. 2001. Two cyclophilin A homologs with shared and distinct functions important for growth and virulence of *Cryptococcus neoformans*. *EMBO Rep.* 2: 511-518.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.

CHROMOSOMAL LOCATION

Genetic locus: PP1B (human) mapping to 15q22.31.

PRODUCT

CyPB siRNA (h2) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CyPB shRNA Plasmid (h2): sc-270696-SH and CyPB shRNA (h2) Lentiviral Particles: sc-270696-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CyPB siRNA (h2) is recommended for the inhibition of CyPB expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

CyPB (k2E2): sc-130626 is recommended as a control antibody for monitoring of CyPB gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CyPB gene expression knockdown using RT-PCR Primer: CyPB (h2)-PR: sc-270696-PR (20 μ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.