

Cdc37 siRNA (h): sc-29255

BACKGROUND

Cell cycle events are regulated by the sequential activation and deactivation of cyclin dependent kinases (Cdks) and by the proteolysis of cyclins. The cell division cycle (Cdc) genes are required at various points in the cell cycle. Cdc25A, Cdc25B and Cdc25C protein Tyrosine phosphatases function as mitotic activators by dephosphorylating Cdc2 p34 on regulatory Tyrosine residues. Cdc6 is the human homolog of *Saccharomyces cerevisiae* Cdc6, which is involved in the initiation of DNA replication. Cdc37 appears to facilitate Cdk4/cyclin D1 complex formation and has been shown to form a stable complex with HSP 90. Cdc34, Cdc27 and Cdc16 function as ubiquitin-conjugating enzymes. Cdc34 is thought to be the structural and functional homolog of *Saccharomyces cerevisiae* Cdc34, which is essential for the G₁ to S phase transition. Cdc16 and Cdc27 are components of the APC (anaphase-promoting complex) which ubiquitinates cyclin B, resulting in cyclin B/Cdk complex degradation.

REFERENCES

- Palmer, R.E., et al. 1990. Mitotic transmission of artificial chromosomes in Cdc mutants of the yeast, *Saccharomyces cerevisiae*. Genetics 125: 763-774.
- Gautier, J., et al. 1991. Cdc25 is a specific Tyrosine phosphatase that directly activates p34^{Cdc2}. Cell 67: 197-211.

CHROMOSOMAL LOCATION

Genetic locus: CDC37 (human) mapping to 19p13.2.

PRODUCT

Cdc37 siRNA (h) is a pool of 4 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Cdc37 shRNA Plasmid (h): sc-29255-SH and Cdc37 shRNA (h) Lentiviral Particles: sc-29255-V as alternate gene silencing products.

For independent verification of Cdc37 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-29255A, sc-29255B, sc-29255C and sc-29255D.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Cdc37 siRNA (h) is recommended for the inhibition of Cdc37 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.


GENE EXPRESSION MONITORING

Cdc37 (E-4): sc-13129 is recommended as a control antibody for monitoring of Cdc37 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).


RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Cdc37 gene expression knockdown using RT-PCR Primer: Cdc37 (h)-PR: sc-29255-PR (20 μ l, 408 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

DATA



Cdc37 siRNA (h): sc-29255. Western blot analysis of Cdc37 expression in non-transfected control (A) and Cdc37 siRNA transfected (B) HeLa cells. Blot probed with Cdc37 (C-11): sc-17758. Actin (I-19): sc-1616 used as specificity and loading control.



Cdc37 siRNA (h): sc-29255. Immunofluorescence staining of methanol-fixed, control HeLa (A) and Cdc37 siRNA silenced HeLa (B) cells showing diminished nuclear staining in the siRNA silenced cells. Cells probed with Cdc37 (H-271): sc-5617.

SELECT PRODUCT CITATIONS

- Babchia, N., et al. 2008. 17-AAG and 17-DMAG-induced inhibition of cell proliferation through B-Raf downregulation in WT B-Raf-expressing uveal melanoma cell lines. Invest. Ophthalmol. Vis. Sci. 49: 2348-2356.
- Kim, H., et al. 2014. Correlation between PDZK1, Cdc37, Akt and breast cancer malignancy: the role of PDZK1 in cell growth through Akt stabilization by increasing and interacting with Cdc37. Mol. Med. 20: 270-279.
- Siddiqui, F.A., et al. 2021. Novel small molecule HSP 90/Cdc37 interface inhibitors indirectly target K-Ras-signaling. Cancers 13: 927.

RESEARCH USE

For research use only, not for use in diagnostic procedures.