

Fyn siRNA (h): sc-29321

BACKGROUND

Src is the human homolog of the v-Src gene of the Rous sarcoma virus, also called avian sarcoma virus or ASV. Src was the first proto-oncogenic non-receptor tyrosine kinase characterized in human. By virtue of common structural motifs, the Src family is composed of nine members in vertebrates, including Src, Yes, Fgr, Frk, Fyn, Lyn, Hck, Lck and Blk. Src-family kinases transduce signals that are involved in the control of a variety of cellular processes, including proliferation, differentiation, motility, and adhesion. Src-family kinases contain an amino terminal cell membrane anchor followed by an SH3 domain and an SH2 domain involved in modular association and activation, respectively. Src-family kinases are normally maintained in an inactive state and can be activated transiently during cellular events such as mitosis. Different subcellular localizations of Src-family kinases may be important for the regulation of specific cellular processes such as mitogenesis, cytoskeletal organization, and membrane trafficking. Fyn and Lck kinases play a key role in T-cell antigen receptor (TCR) signaling. The human Fyn gene maps to chromosome 6q21 and encodes a 537 amino acid protein.

CHROMOSOMAL LOCATION

Genetic locus: FYN (human) mapping to 6q21.

PRODUCT

Fyn siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Fyn shRNA Plasmid (h): sc-29321-SH and Fyn shRNA (h) Lentiviral Particles: sc-29321-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Fyn siRNA (h) is recommended for the inhibition of Fyn expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Fyn (15): sc-434 is recommended as a control antibody for monitoring of Fyn gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Fyn gene expression knockdown using RT-PCR Primer: Fyn (h)-PR: sc-29321-PR (20 μ l, 542 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Qiao, F., et al. 2008. Distinct functions of Cdk5(Y15) phosphorylation and Cdk5 activity in stress fiber formation and organization. *Exp. Cell Res.* 314: 3542-3550.
- Huang, Y.T., et al. 2008. Sphingosine 1-phosphate induces platelet/endothelial cell adhesion molecule-1 phosphorylation in human endothelial cells through c-Src and Fyn. *Cell. Signal.* 20: 1521-1527.
- Sinha, S., et al. 2009. Dopamine regulates phosphorylation of VEGF receptor 2 by engaging Src-homology-2-domain-containing protein tyrosine phosphatase 2. *J. Cell Sci.* 122: 3385-3392.
- Knezevic, N., et al. 2009. The G protein β subunit mediates reannealing of adherens junctions to reverse endothelial permeability increase by thrombin. *J. Exp. Med.* 206: 2761-2777.
- Xue, Q., et al. 2010. *Bacillus anthracis* spore entry into epithelial cells is an Actin-dependent process requiring c-Src and PI3K. *PLoS ONE* 5: e11665.
- Yadav, V. and Denning, M.F. 2011. Fyn is induced by Ras/PI3K/Akt signaling and is required for enhanced invasion/migration. *Mol. Carcinog.* 50: 346-352.
- Liu, W., et al. 2012. LKB1/STK11 inactivation leads to expansion of a prometastatic tumor subpopulation in melanoma. *Cancer Cell* 21: 751-764.
- Kim, A.Y., et al. 2012. Enhanced antioxidant effect of prenylated polypeptols as Fyn inhibitor. *Free Radic. Biol. Med.* 53: 1198-1208.
- Ninio-Many, L., et al. 2013. MicroRNA-125a-3p reduces cell proliferation and migration by targeting Fyn. *J. Cell Sci.* 126: 2867-2876.
- Grossman, H., et al. 2015. A novel regulatory pathway in granulosa cells, the LH/human chorionic gonadotropin-microRNA-125a-3p-Fyn pathway, is required for ovulation. *FASEB J.* 29: 3206-3216.
- Xie, Y.G., et al. 2016. FYN promotes breast cancer progression through epithelial-mesenchymal transition. *Oncol. Rep.* 36: 1000-1006.
- Venkata Subbaiah, K.C., et al. 2019. hnRNP L-mediated RNA switches function as a hypoxia-induced translational regulon. *Biochem. Biophys. Res. Commun.* 516: 753-759.

RESEARCH USE

For research use only, not for use in diagnostic procedures.