

Lyn siRNA (h): sc-29393

BACKGROUND

Src is the human homolog of the v-Src gene of the Rous sarcoma virus, also known as avian sarcoma virus or ASV. Src is the first proto-oncogenic non-receptor tyrosine kinase characterized in human. By virtue of common structural motifs, the Src family is composed of nine members in vertebrates, including Src, Yes, Fgr, Frk, Fyn, Lyn, Hck, Lck and Blk. Src-family kinases transduce signals that are involved in the control of a variety of cellular processes, including proliferation, differentiation, motility and adhesion. Src-family kinases contain an amino-terminal cell membrane anchor followed by an SH3 domain and an SH2 domain involved in modular association and activation, respectively. Src-family kinases are normally maintained in an inactive state and can be activated transiently during cellular events such as mitosis. The human Lyn gene maps to chromosome 8q12.1 and encodes a 505 amino acid protein. Hematopoietic tissues predominantly express Lyn, which influences normal immunoglobulin production and regulation.

CHROMOSOMAL LOCATION

Genetic locus: LYN (human) mapping to 8q12.1.

PRODUCT

Lyn siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Lyn shRNA Plasmid (h): sc-29393-SH and Lyn shRNA (h) Lentiviral Particles: sc-29393-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Lyn siRNA (h) is recommended for the inhibition of Lyn expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Lyn (H-6): sc-7274 is recommended as a control antibody for monitoring of Lyn gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Lyn gene expression knockdown using RT-PCR Primer: Lyn (h)-PR: sc-29393-PR (20 μ l, 548 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

1. Basu, S., et al. 2008. *Helicobacter pylori* protein HP0175 transactivates epidermal growth factor receptor through TLR4 in gastric epithelial cells. *J. Biol. Chem.* 283: 32369-32376.
2. Kim, M.J., et al. 2008. Activation of Lck is critically required for sphingosine-induced conformational activation of Bak and mitochondrial cell death. *Biochem. Biophys. Res. Commun.* 370: 353-358.
3. Su, N., et al. 2012. Lyn is involved in CD24-induced ERK1/2 activation in colorectal cancer. *Mol. Cancer* 11: 43.
4. Li, G., et al. 2013. Lyn mitigates mouse airway remodeling by downregulating the TGF- β 3 isoform in house dust mite models. *J. Immunol.* 191: 5359-5370.
5. Zardan, A., et al. 2014. Lyn tyrosine kinase regulates androgen receptor expression and activity in castrate-resistant prostate cancer. *Oncogenesis* 3: e115.
6. Chen, Y.J., et al. 2017. Hydroquinone-induced FOXp3-ADAM17-Lyn-Akt-p21 signaling axis promotes malignant progression of human leukemia U937 cells. *Arch. Toxicol.* 91: 983-997.
7. Wang, X., et al. 2017. Lyn regulates mucus secretion and MUC5AC via the STAT6 signaling pathway during allergic airway inflammation. *Sci. Rep.* 7: 42675.
8. Wang, X., et al. 2017. Lyn kinase represses mucus hypersecretion by regulating IL-13-induced endoplasmic reticulum stress in asthma. *EBioMedicine* 15: 137-149.
9. Liang, X., et al. 2019. Lyn regulates epithelial-mesenchymal transition in CS-exposed model through Smad2/3 signaling. *Respir. Res.* 20: 201.
10. Yuan, M., et al. 2020. DC-SIGN-LEF1/TCF1-miR-185 feedback loop promotes colorectal cancer invasion and metastasis. *Cell Death Differ.* 27: 379-395.
11. Seok, H.J., et al. 2021. Novel miR-5088-5p promotes malignancy of breast cancer by inhibiting DBC2. *Mol. Ther. Nucleic Acids* 25: 127-142.

RESEARCH USE

For research use only, not for use in diagnostic procedures.