

# nm23-H1 siRNA (m): sc-29415

## BACKGROUND

The nm23 gene, a potential suppressor of metastasis, was originally identified by differential hybridization between two murine melanoma sub-lines, one with a high and the second with a low metastatic capacity. Highly metastatic sub-lines exhibit much lower levels of nm23 than less metastatic cells. Based on sequence analysis, nm23 appears highly related to nucleotide diphosphate kinases (NDP). In humans, NDP kinases A and B are identical to two isoforms of human nm23 homologs, namely nm23-H1 and H2, respectively. nm23-H2 is identical in sequence to PuF, a transcription factor that binds to nuclease-hypersensitive elements at positions 142 to 115 of the human c-Myc promoter.

## REFERENCES

1. Steeg, P.S., et al. 1988. Evidence for a novel gene associated with low tumor metastatic potential. *J. Natl. Cancer Inst.* 80: 200-209.
2. Lacombe, M., et al. 1990. Functional cloning of a nucleoside diphosphate kinase from *Dictyostelium discoideum*. *J. Biol. Chem.* 265: 10012-10018.
3. Kimura, N., et al. 1990. Isolation and characterization of a cDNA clone encoding rat nucleoside diphosphate kinase. *J. Biol. Chem.* 265: 15744-15749.
4. Stahl, J.A., et al. 1991. Identification of a second human nm23 gene, nm23-H2. *Cancer Res.* 51: 445-449.
5. Urano, T., et al. 1992. Molecular cloning and functional expression of the second mouse nm23/NDP kinase gene, nm23-M2. *FEBS Lett.* 309: 358-362.
6. Urano, T., et al. 1993. Expression of nm23/NDP kinase proteins on the cell surface. *Oncogene* 8: 1371-1376.
7. Postel, E.H., et al. 1993. Human c-Myc transcription factor PuF identified as nm23-H2 nucleoside diphosphate kinase, a candidate suppressor of tumor metastasis. *Science* 261: 478-480.

## CHROMOSOMAL LOCATION

Genetic locus: Nme1 (mouse) mapping to 11 D.

## PRODUCT

nm23-H1 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see nm23-H1 shRNA Plasmid (m): sc-29415-SH and nm23-H1 shRNA (m) Lentiviral Particles: sc-29415-V as alternate gene silencing products.

For independent verification of nm23-H1 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-29415A, sc-29415B and sc-29415C.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

nm23-H1 siRNA (m) is recommended for the inhibition of nm23-H1 expression in mouse cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

nm23-H1 (NM301): sc-465 is recommended as a control antibody for monitoring of nm23-H1 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor nm23-H1 gene expression knockdown using RT-PCR Primer: nm23-H1 (m)-PR: sc-29415-PR (20  $\mu$ l). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.