

Smad4 siRNA (h): sc-29484

BACKGROUND

Smad proteins, the mammalian homologs of the *Drosophila* mothers against dpp (Mad) have been implicated as downstream effectors of TGF β /BMP signaling. Smad1 (also designated Madr1 or JV4-1), Smad5 and mammalian Smad8 (also designated Smad9 or MADH6) are effectors of BMP2 and BMP4 function while Smad2 (also designated Madr2 or JV18-1) and Smad3 are involved in TGF β and Activin-mediated growth modulation. Smad4 (also designated DPC4) has been shown to mediate all of the above activities through interaction with various Smad family members. Smad6 and Smad7 regulate the response to Activin/TGF β signaling by interfering with TGF β -mediated phosphorylation of other Smad family members.

CHROMOSOMAL LOCATION

Genetic locus: SMAD4 (human) mapping to 18q21.2.

PRODUCT

Smad4 siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Smad4 shRNA Plasmid (h): sc-29484-SH and Smad4 shRNA (h) Lentiviral Particles: sc-29484-V as alternate gene silencing products.

For independent verification of Smad4 (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-29484A, sc-29484B and sc-29484C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Smad4 siRNA (h) is recommended for the inhibition of Smad4 expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

Smad4 (B-8): sc-7966 is recommended as a control antibody for monitoring of Smad4 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Smad4 gene expression knockdown using RT-PCR Primer: Smad4 (h)-PR: sc-29484-PR (20 μ l, 555 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

1. Kennard, S., et al. 2008. Transforming growth factor- β (TGF-1) down-regulates Notch3 in fibroblasts to promote smooth muscle gene expression. *J. Biol. Chem.* 283: 1324-1333.
2. Ren, S., et al. 2009. Transforming growth factor- β 2 upregulates sphingosine kinase-1 activity, which in turn attenuates the fibrotic response to TGF- β 2 by impeding CTGF expression. *Kidney Int.* 76: 857-867.
3. Pan, H., et al. 2010. Smad4 is required for development of maximal endotoxin tolerance. *J. Immunol.* 184: 5502-5509.
4. Wang, H., et al. 2014. Hepcidin is directly regulated by Insulin and plays an important role in iron overload in streptozotocin-induced diabetic rats. *Diabetes* 63: 1506-1518.
5. Chatterjee, A., et al. 2015. Rapamycin-induced G₁ cell cycle arrest employs both TGF- β and Rb pathways. *Cancer Lett.* 360: 134-140.
6. Kim, R.Y., et al. 2015. Estrogen modulates bone morphogenetic protein-induced sclerostin expression through the Wnt signaling pathway. *Tissue Eng. Part A* 21: 2076-2088.
7. Hossain, D.M., et al. 2015. MEK inhibition prevents tumour-shed transforming growth factor- β -induced T-regulatory cell augmentation in tumour milieu. *Immunology* 144: 561-573.
8. Wa, Q., et al. 2018. Downregulation of miR-19a-3p promotes invasion, migration and bone metastasis via activating TGF- β signaling in prostate cancer. *Oncol. Rep.* 39: 81-90.
9. Yang, X.H., et al. 2018. MiR-3147 serves as an oncomiR in vulvar squamous cell cancer via Smad4 suppression. *Mol. Med. Rep.* 17: 6397-6404.
10. Lou, Z., et al. 2018. ESE-1 suppresses the growth, invasion and migration of human NSCLC cells and tumor formation *in vivo*. *Oncol. Rep.* 40: 1734-1742.
11. Chao, D., et al. 2019. AZD3759 induces apoptosis in hepatoma cells by activating a p53-Smad4 positive feedback loop. *Biochem. Biophys. Res. Commun.* 509: 535-540.

RESEARCH USE

For research use only, not for use in diagnostic procedures.