

Ah Receptor siRNA (h): sc-29654

BACKGROUND

2,3,7,8-tetrachlorodibenzo-p-dioxin (TCDD) is the prototype for a family of toxic halogenated aromatic compounds that are thought to cause adverse reproductive, immunologic and metabolic effects. Many biological responses to TCDD are mediated through ligand binding to the aromatic hydrocarbon (Ah) receptor, also known as AhR. Ah Receptor is a ligand dependent transcription factor that interacts with specific DNA sequences, termed xenobiotic responsive elements (XREs), and that lies upstream of TCDD-inducible genes. Upon binding to the ligand, Ah Receptor binds to the Ah Receptor nuclear translocator (Arnt), and the complex is translocated from the cytoplasm to the nucleus. Arnt is required for Ah Receptor to bind to XRE. Ah Receptor and Arnt are members of a family of transcription factors that contain a basic helix-loop-helix motif and a common "PAS" motif.

CHROMOSOMAL LOCATION

Genetic locus: AHR (human) mapping to 7p21.1.

PRODUCT

Ah Receptor siRNA (h) is a target-specific 19-25 nt siRNA designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see Ah Receptor shRNA Plasmid (h): sc-29654-SH and Ah Receptor shRNA (h) Lentiviral Particles: sc-29654-V as alternate gene silencing products.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

Ah Receptor siRNA (h) is recommended for the inhibition of Ah Receptor expression in human cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

GENE EXPRESSION MONITORING

Ah Receptor (A-3): sc-133088 is recommended as a control antibody for monitoring of Ah Receptor gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor Ah Receptor gene expression knockdown using RT-PCR Primer: Ah Receptor (h)-PR: sc-29654-PR (20 μ l, 442 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Kasao, A., et al. 2006. High levels of Dioxin-like potential in cigarette smoke evidenced by *in vitro* and *in vivo* biosensing. *Cancer Res.* 66: 7143-7150.
- Yisireyili, M., et al. 2014. Indoxyl sulfate-induced activation of (pro)renin receptor promotes cell proliferation and tissue factor expression in vascular smooth muscle cells. *PLoS ONE* 9: e109268.
- Kumar, R., et al. 2015. Quercetin-6-C- β -D-glucopyranoside, natural analog of quercetin exhibits anti-prostate cancer activity by inhibiting Akt-mTOR pathway via aryl hydrocarbon receptor. *Biochimie* 119: 68-79.
- Ng, H.Y., et al. 2016. Indoxyl sulfate downregulates mas receptor via aryl hydrocarbon receptor/nuclear factor- κ B, and induces cell proliferation and tissue factor expression in vascular smooth muscle cells. *Nephron* 133: 205-212.
- Gao, W.Q., et al. 2017. Paclitaxel-mediated human aryl hydrocarbon receptor mRNA translation by an internal ribosomal entry site-dependent mechanism. *Oncol. Rep.* 38: 3211-3219.
- Bogoevska, V., et al. 2017. HRG/HER2/HER3 signaling promotes AhR-mediated Memo-1 expression and migration in colorectal cancer. *Oncogene* 36: 2394-2404.
- Das, D.N., et al. 2017. DNA damage by 2,3,7,8-tetrachlorodibenzo-p-Dioxin-induced p53-mediated apoptosis through activation of cytochrome P450/aryl hydrocarbon receptor. *Environ. Toxicol. Pharmacol.* 55: 175-185.
- Solis, N.V., et al. 2017. The aryl hydrocarbon receptor governs epithelial cell invasion during oropharyngeal candidiasis. *mBio* 8: e00025-17.
- Yabuuchi, N., et al. 2021. Suppressed hepatic production of indoxyl sulfate attenuates cisplatin-induced acute kidney injury in sulfotransferase 1a1-deficient mice. *Int. J. Mol. Sci.* 22: 1764.
- Zhang, X., et al. 2021. Blockade of IDO-kynurenine-AhR axis ameliorated colitis-associated colon cancer via inhibiting immune tolerance. *Cell. Mol. Gastroenterol. Hepatol.* E-published.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.