



ATF-3 siRNA (m): sc-29758

BACKGROUND

Eukaryotic gene transcription is regulated by sequence-specific transcription factors which bind modular *cis* acting promotor and enhancer elements. The cAMP response element (CRE), one of the best studied of such elements, consists of the palindromic octanucleotide TGACGTCA. Several CRE binding proteins have been identified within the ATF/CREB family, the best characterized of which include CREB-1, CREB-2 (also designated ATF-4), ATF-1, ATF-2 and ATF-3. These proteins share highly related COOH-terminal leucine zipper dimerization and basic DNA binding domains but are highly divergent in their amino terminal domains. Although each of the ATF/CREB proteins appear capable of binding CRE in its homodimeric form, certain of these also bind as heterodimers, both within the ATF/CREB family and even with members of the AP-1 transcription factor family.

CHROMOSOMAL LOCATION

Genetic locus: Atf3 (mouse) mapping to 1 H6.

PRODUCT

ATF-3 siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see ATF-3 shRNA Plasmid (m): sc-29758-SH and ATF-3 shRNA (m) Lentiviral Particles: sc-29758-V as alternate gene silencing products.

For independent verification of ATF-3 (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-29758A, sc-29758B and sc-29758C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

ATF-3 siRNA (m) is recommended for the inhibition of ATF-3 expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

GENE EXPRESSION MONITORING

ATF-3 (A-8): sc-518032 is recommended as a control antibody for monitoring of ATF-3 gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor ATF-3 gene expression knockdown using RT-PCR Primer: ATF-3 (m)-PR: sc-29758-PR (20 μ l, 554 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

SELECT PRODUCT CITATIONS

- Kim, J.Y., et al. 2010. Chronic ethanol consumption-induced pancreatic β -cell dysfunction and apoptosis through glucokinase nitration and its down-regulation. *J. Biol. Chem.* 285: 37251-37262.
- Jang, M.K., et al. 2011. ATF3 represses PDX-1 expression in pancreatic β -cells. *Biochem. Biophys. Res. Commun.* 412: 385-390.
- Kim, J., et al. 2014. Metformin suppresses lipopolysaccharide (LPS)-induced inflammatory response in murine macrophages via activating transcription factor-3 (ATF-3) induction. *J. Biol. Chem.* 289: 23246-23255.
- Bae, Y.A. and Cheon, H.G. 2016. Activating transcription factor-3 induction is involved in the anti-inflammatory action of berberine in RAW 264.7 murine macrophages. *Korean J. Physiol. Pharmacol.* 20: 415-424.
- Qian, L., et al. 2017. Activating transcription factor 3 (ATF-3) protects against lipopolysaccharide-induced acute lung injury via inhibiting the expression of TL1A. *J. Cell. Physiol.* 232: 3727-3734.
- Hu, C., et al. 2017. Frontline science: ATF3 is responsible for the inhibition of TNF- α release and the impaired migration of acute ethanol-exposed monocytes and macrophages. *J. Leukoc. Biol.* 101: 633-642.
- Hernandez, G., et al. 2020. Pancreatitis is an FGF21-deficient state that is corrected by replacement therapy. *Sci. Transl. Med.* 12: eaay5186.
- An, M.Y., et al. 2020. Anti-inflammatory effects of leaf and branch extracts of honeyberry (*Lonicera caerulea*) on lipopolysaccharide-stimulated RAW 264.7 cells through ATF3 and Nrf2/HO-1 activation. *Mol. Med. Rep.* 22: 5219-5230.
- Saha, S., et al. 2021. *Leishmania donovani* targets host transcription factor NRF2 to activate anti-oxidant enzyme HO-1 and transcriptional repressor ATF3 for establishing infection. *Infect. Immun.* 89: e0076420.
- Aguilar-Recarte, D., et al. 2022. A positive feedback loop between AMPK and GDF15 promotes metformin antidiabetic effects. *Pharmacol. Res.* 187: 106578.
- Zhang, Y., et al. 2024. MSC-derived mitochondria promote axonal regeneration via Atf3 gene up-regulation by ROS induced DNA double strand breaks at transcription initiation region. *Cell Commun. Signal.* 22: 240.

RESEARCH USE

For research use only, not for use in diagnostic procedures.