



# CDP siRNA (h): sc-35051

## BACKGROUND

CDP (for CCAAT displacement protein) has been identified as a repressor for transcription of developmentally regulated genes. It is a homeodomain protein that appears to compete with transcriptional activating proteins for binding to the promoter regions of various genes. CDP contains three cut repeats which function as DNA binding domains. It has been demonstrated that cut repeat domains have the capacity to bind to DNA in conjunction with or independently of homeodomain DNA binding. CDP has been shown to be the DNA-binding subunit of the HiNF-D complex, which contains cyclin A, Cdc2 and an Rb-related protein, in addition to CDP. Histone expression is required for the transition to S phase in the cell cycle. The HiNF-D complex regulates the transcription of histone H4, H3 and H1 genes, allowing cells to progress from G<sub>1</sub> to S phase.

## REFERENCES

1. Stein, G.S., et al, eds. 1984. Histone genes. New York: Wiley.
2. Neufeld, E.J., et al. 1992. Human CCAAT displacement protein is homologous to the *Drosophila* homeoprotein, CUT. Nat. Genet. 1: 50-55.
3. Valarche, I., et al. 1993. The mouse homeodomain protein Phox2 regulates NCAM promoter activity in concert with Cux/CDP and is a putative determinant of neurotransmitter phenotype. Development 119: 881-896.
4. Harada, R., et al. 1994. Conserved CUT repeats in the human CUT homeodomain protein function as DNA binding domains. J. Biol. Chem. 269: 2062-2067.

## CHROMOSOMAL LOCATION

Genetic locus: CUX1 (human) mapping to 7q22.1.

## PRODUCT

CDP siRNA (h) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10  $\mu$ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CDP shRNA Plasmid (h): sc-35051-SH and CDP shRNA (h) Lentiviral Particles: sc-35051-V as alternate gene silencing products.

For independent verification of CDP (h) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-35051A, sc-35051B and sc-35051C.

## STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNases and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330  $\mu$ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330  $\mu$ l of RNase-free water makes a 10  $\mu$ M solution in a 10  $\mu$ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

## APPLICATIONS

CDP siRNA (h) is recommended for the inhibition of CDP expression in human cells.

## SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10  $\mu$ M in 66  $\mu$ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

## GENE EXPRESSION MONITORING

CDP (B-10): sc-514008 is recommended as a control antibody for monitoring of CDP gene expression knockdown by Western Blotting (starting dilution 1:200, dilution range 1:100-1:1000) or immunofluorescence (starting dilution 1:50, dilution range 1:50-1:500).

To ensure optimal results, the following support reagents are recommended: 1) Western Blotting: use m-IgG $\kappa$  BP-HRP: sc-516102 or m-IgG $\kappa$  BP-HRP (Cruz Marker): sc-516102-CM (dilution range: 1:1000-1:10000), Cruz Marker™ Molecular Weight Standards: sc-2035, UltraCruz® Blocking Reagent: sc-516214 and Western Blotting Luminol Reagent: sc-2048. 2) Immunofluorescence: use m-IgG $\kappa$  BP-FITC: sc-516140 or m-IgG $\kappa$  BP-PE: sc-516141 (dilution range: 1:50-1:200) with UltraCruz® Mounting Medium: sc-24941 or UltraCruz® Hard-set Mounting Medium: sc-359850.

## RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CDP gene expression knockdown using RT-PCR Primer: CDP (h)-PR: sc-35051-PR (20  $\mu$ l, 405 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

## SELECT PRODUCT CITATIONS

1. Li, J., et al. 2007. Protein kinase C-mediated modulation of FIH-1 expression by the homeodomain protein CDP/Cut/Cux. Mol. Cell. Biol. 27: 7345-7353.
2. Wang, J., et al. 2017. CUTL1 induces epithelial-mesenchymal transition in non-small cell lung cancer. Oncol. Rep. 37: 3068-3074.

## RESEARCH USE

For research use only, not for use in diagnostic procedures.

## PROTOCOLS

See our web site at [www.scbt.com](http://www.scbt.com) for detailed protocols and support products.