



CRABP-I siRNA (m): sc-35104

BACKGROUND

The cellular retinoic acid-binding protein (CRABP)-I and a related isoform, CRABP-II, bind retinoic acid (RA), an important regulator of cell growth and differentiation in fetal and adult tissues. These CRABP proteins mediate the downstream effects of RA in distinct ways. CRABP-I negatively regulates the activity of RA by enhancing the production of RA-metabolizing enzymes and increasing the rate at which RA is degraded. CRABP-II enhances the effects of RA by directly interacting with RA receptors (RAR) and, in turn, promoting the formation of RAR-RA complexes and stimulating RA-mediated gene transcription. Both CRABP-I and CRABP-II are expressed in the embryo, and CRABP-I is ubiquitously expressed in various adult tissues. The expression of CRABP-II is elevated in cells that synthesize relatively large amounts of RA, and it is also predominantly expressed in skin, uterus, ovary and in the choroid plexus.

REFERENCES

1. Wei, L.N., Tsao, J.L., Chu, Y.S., Jeannotte, L. and Nguyen-Huu, M.C. 1990. Molecular cloning and transcriptional mapping of the mouse cellular retinoic acid-binding protein gene. *DNA Cell Biol.* 9: 471-478.
2. Giguere, V., Lyn, S., Yip, P., Siu, C.H. and Amin, S. 1990. Molecular cloning of cDNA encoding a second cellular retinoic acid-binding protein. *Proc. Natl. Acad. Sci. USA* 87: 6233-6237.
3. Boylan, J.F. and Gudas, L.J. 1992. The level of CRABP-I expression influences the amounts and types of all-*trans*-retinoic acid metabolites in F9 teratocarcinoma stem cells. *J. Biol. Chem.* 267: 21486-21491.
4. Gorry, P., Lufkin, T., Dierich, A., Rochette-Egly, C., Decimo, D., Dolle, P., Mark, M., Durand, B. and Chambon, P. 1994. The cellular retinoic acid binding protein I is dispensable. *Proc. Natl. Acad. Sci. USA* 91: 9032-9036.
5. Astrom, A., Pettersson, U., Chambon, P. and Voorhees, J.J. 1994. Retinoic acid induction of human cellular retinoic acid-binding protein-II gene transcription is mediated by retinoic acid receptor-retinoid X receptor heterodimers bound to one far upstream retinoic acid-responsive element with 5-base pair spacing. *J. Biol. Chem.* 269: 22334-22339.
6. Zheng, W.L., Bucco, R.A., Schmitt, M.C., Wardlaw, S.A., and Ong, D.E. 1996. Localization of cellular retinoic acid-binding protein (CRABP) II and CRABP in developing rat testis. *Endocrinology* 137: 5028-5035.

CHROMOSOMAL LOCATION

Genetic locus: Crabp1 (mouse) mapping to 9 A5.3.

PRODUCT

CRABP-I siRNA (m) is a pool of 3 target-specific 19-25 nt siRNAs designed to knock down gene expression. Each vial contains 3.3 nmol of lyophilized siRNA, sufficient for a 10 μ M solution once resuspended using protocol below. Suitable for 50-100 transfections. Also see CRABP-I shRNA Plasmid (m): sc-35104-SH and CRABP-I shRNA (m) Lentiviral Particles: sc-35104-V as alternate gene silencing products.

For independent verification of CRABP-I (m) gene silencing results, we also provide the individual siRNA duplex components. Each is available as 3.3 nmol of lyophilized siRNA. These include: sc-35104A, sc-35104B and sc-35104C.

STORAGE AND RESUSPENSION

Store lyophilized siRNA duplex at -20° C with desiccant. Stable for at least one year from the date of shipment. Once resuspended, store at -20° C, avoid contact with RNAses and repeated freeze thaw cycles.

Resuspend lyophilized siRNA duplex in 330 μ l of the RNase-free water provided. Resuspension of the siRNA duplex in 330 μ l of RNase-free water makes a 10 μ M solution in a 10 μ M Tris-HCl, pH 8.0, 20 mM NaCl, 1 mM EDTA buffered solution.

APPLICATIONS

CRABP-I siRNA (m) is recommended for the inhibition of CRABP-I expression in mouse cells.

SUPPORT REAGENTS

For optimal siRNA transfection efficiency, Santa Cruz Biotechnology's siRNA Transfection Reagent: sc-29528 (0.3 ml), siRNA Transfection Medium: sc-36868 (20 ml) and siRNA Dilution Buffer: sc-29527 (1.5 ml) are recommended. Control siRNAs or Fluorescein Conjugated Control siRNAs are available as 10 μ M in 66 μ l. Each contain a scrambled sequence that will not lead to the specific degradation of any known cellular mRNA. Fluorescein Conjugated Control siRNAs include: sc-36869, sc-44239, sc-44240 and sc-44241. Control siRNAs include: sc-37007, sc-44230, sc-44231, sc-44232, sc-44233, sc-44234, sc-44235, sc-44236, sc-44237 and sc-44238.

RT-PCR REAGENTS

Semi-quantitative RT-PCR may be performed to monitor CRABP-I gene expression knockdown using RT-PCR Primer: CRABP-I (m)-PR: sc-35104-PR (20 μ l, 437 bp). Annealing temperature for the primers should be 55-60° C and the extension temperature should be 68-72° C.

RESEARCH USE

For research use only, not for use in diagnostic procedures.

PROTOCOLS

See our web site at www.scbt.com for detailed protocols and support products.